Instructions for Use Cannula Insertion System

The Cannula Insertion System is intended for use to establish an entry site conduit for passing ophthalmic instruments and fluids used to perform posterior ophthalmic surgical procedures.

REF	Gauge
5700GEN / 5700GST	23GA
5600GEN / 5600GST	25GA
5500GEN / 5500GST	27GA

Instructions for Use

1. Fixate the conjunctiva in place during the insertion procedure.
2. Insert the cannula inserter into the sclera.
3. Remove the blade and fix the cannula on the eye.
4. Retract the blade.
5. Connect the infusion tubing into the valve and switch irrigation on.
6. Repeat step. 1 to 4 for the required number of ports.
7. Remove the cannulas when surgery is complete.

INDICATION

The Cannula Insertion System is indicated for use in most vitreoretinal surgical procedures which utilize a pars plana entry site.

CONTRAINDICATION

Do not use in the presence of poorly adhesive intraocular tumors such as retinoblastomas. Alternative methods should be considered in case of: severe scleral thinning; wounds from recent ocular surgery; pediatric anatomy.

⚠ WARNING

- Risk of cross-contamination if product is reused.
- Do not use if blade is damaged.
- Do not actuate cutting instrument while passing through the valve.
- Ensure proper placement of trocar cannulas to prevent sub-retinal infusion.
- Leaking sclerotomy may lead to post-operative hypotony.
- This device is for short-term use during surgical procedures. Do not use for more than 24 hours.

⚠ CAUTION

Federal (USA) law restricts this device to use or sale by or on order of a physician.

INTENDED USERS

The Cannula Insertion System is intended to be used by healthcare professionals.

EN

STERILE EO Rx Only

HOYA Surgical Optics GmbH
10 De-Saint-Exupéry Strasse
60549 Frankfurt
Germany
EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

TARGET GROUP

Because of the wide variety of potential indications, the patient population is limited only by the physician judgement of suitability for surgery. Populations with known limitations on cannula usage include infants (particularly less than four years of age) and highly myopic patients.

CLINICAL BENEFIT

The Cannula insertion System reduces surgical time, reduces post-surgical discomfort and reduces surgically-induced astigmatism.

LIMITATION

Only use with instruments that are compatible in size.

If any MID Labs products are received in a defective condition, your distributor or MID Labs is to be notified immediately.

Medical Instrument Development
Laboratories, Inc.
557 McCormick Street,
San Leandro, CA 94577 U.S.A.

Phone: 1-800-929-5227 (Inside U.S. only)
Phone: 510-357-3952
Fax: 510-357-1582
Email: info@midlabs.com

Each MID Labs product is identified by a lot number. Refer to this number when discussing the product with your MID Labs representative. Any serious incident that has occurred in relation to this device should be reported to MID Labs and the competent authority of the Member State in which the user and/or patient is established.

Glossary of Symbols			
Symbol	Title/Description	Applicable Standard	Reference Number
	Manufacturer	EN ISO 15223-1	5.1.1
	Authorized Representative in the European Community	EN ISO 15223-1	5.1.2
	Date of Manufacture	EN ISO 15223-1	5.1.3
	Use-by date	EN ISO 15223-1	5.1.4
	Batch code	EN ISO 15223-1	5.1.5
	Catalogue number	EN ISO 15223-1	5.1.6
	Importer	EN ISO 15223-1	5.1.8
	Sterilized using ethylene oxide	EN ISO 15223-1	5.2.3
	Do not re-sterilize	EN ISO 15223-1	5.2.6
	Do not use if package is damaged	EN ISO 15223-1	5.2.8
	Single sterile barrier system	EN ISO 15223-1	5.2.11
	Single sterile barrier system with protective packaging outside	EN ISO 15223-1	5.2.14
	Do not re-use	EN ISO 15223-1	5.4.2
	Caution	EN ISO 15223-1	5.4.4
	Medical Device	EN ISO 15223-1	5.7.7
	Translation	EN ISO 15223-1	5.7.8
	Consult Instructions for use	EN ISO 15223-1	5.4.3
	Caution: Federal (U.S.) Law Restricts this device to sale by or on the order of a physician.	21 CFR 801	§ 801.109 (b)(1)
	European Conformity - CE Mark	Council Directive 93/42/EEC	ANNEX XII

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Medical Instrument Development Laboratories, Inc.