Directions for Use
Bi-Blade® Vitreous Cutter With Tubing Set

The MID Labs Bi-Blade® Vitreous Cutter with Tubing Set is intended for use to remove vitreous and intra-ocular tissue, and may be used with any compatible ophthalmic surgical system producing minimum actuation and maximum residual pressure at tubing end as specified in the system requirements table below.

System Requirements:

Table with 4 columns: REF, Minimum Actuation Pressure, Maximum Residual Pressure, MID Labs Systems. Row 1: 28- -BB, 269KPa (39.0 psi), 152 KPa (22.0 psi), ---

STERILE R icons: No re-sterilize, No re-use, Do not use if package is damaged, PHT DEHP

CAUTION: Do not operate above 2500 cycles per minute equivalent to 5000 cuts per minute with Bi-Blade® technology.

CAUTION: Bi-Blade® aspiration flow may be greater than expected. Start with a reduced vacuum level.

WARNING: Sterility not guaranteed if package has been opened or damaged. This product is for single patient use only. Do not re-sterilize. Risk of cross-contamination if product is re-used.

CAUTION: Federal (USA) law restricts this device to use or sale by or on order of a physician.

CAUTION: 344 KPa (50 psi) Maximum Pressure at Tubing End


WARNING: Test Vitreous Cutter before use. If bubbles are seen at cutting port, do not use cutter.

CAUTION: Do not operate the cutter dry since the cutting action is dependent on a liquid environment.

CAUTION: Failure to properly attach the air supply tubing to the appropriate pressure source will affect the operation of the vitrectomy cutter.

INSTALLATION

- 1. Tubing set is color coded for easy identification; GREEN striped line for cutting actuation pressure and CLEAR tubing for aspiration.
- 2. Handle using appropriate aseptic technique.
- 3. The GREEN striped tubing is to be attached to the vitrectomy actuation port on the system console. Disconnect the plug from the tubing set and attach it to the GREY actuation port with the icon on the system console. Attach the connector with GREEN striped tubing to the BLACK actuation port with the icon on the system console.
- 4. The CLEAR tubing is to be attached to the aspiration system of vitrectomy system console.

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Table with 4 columns: Símbolo, Título/descripción, Norma correspondiente, Número de referencia. Rows include Manufacturer, Authorized Representative, Use-by date, Batch code, Catalogue number, Sterilized using irradiation, Do not re-sterilize, Do not re-use, Do not use if package is damaged, Consult instructions for use, Caution, Rx Only, European Conformity, Contains or presence of phthalates.

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Medical Instrument Development Laboratories, Inc.
San Leandro, CA 94577 U.S.A.
Tel: 1-800-929-5227 (solo U.S.A.)