INSTRUCTIONS FOR USE

BioComFold

FOLDABLE INTRAOCULAR LENSES

Manufacturer

MORCHER® GmbH

Kapuzinerweg 12

70374 Stuttgart

GERMANY

CE0483

F41042V2021-08

GRAPHICS 4 - 5

ENGLISH

INSTRUCTIONS FOR USE 6 - 15

DEUTSCH

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ESPAÑOL




INSTRUCCIONES DE USO..... 26 - 35




FRANÇAIS

NOTICE..... 36 - 45

ITALIANO

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INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

The MORCHER® BioComFold (BCF) intraocular lens (IOL), hereafter referred to as lens or implant, is a sterile optical medical implant designed to replace or supplement the natural lens. The quality of the material used permits folding the optic for small incision insertion. The lenses consist of a hydrophilic material (acrylic copolymer) including a copolymerized UV filter and, if applicable, an additional copolymerized HEV filter (blue light filter). The lens is available with different size optics. Each lens is tested separately for its optical quality and refractive power.

The BioComFold IOLs are available with different diameters and in various forms in order to ensure optimal fixing of the implant in the eye.

Accessory: Caliper Ring

The Caliper Ring is an accessory designed to size and centre the position of the anterior capsulorhexis. The rings are available in different sizes.

INTENDED USE

The MORCHER® BioComFold (BCF) intraocular lens is intended to serve as replacement of or complement to natural lens. The lenses are available in different modifications:

- For fixation into the capsular bag

They are used after ECCE and ICCE (extra-capsular or intra-capsular cataract extraction).

INDICATIONS

- For senile cataracts
- As correction of a phakic or aphakic eye where other optical aids are not available or are not implementable
- For traumatic cataracts
- In exceptional situations with congenital or early-acquired cataracts

CONTRAINDICATIONS

The usual contraindications for a cataract operation apply here, as well as the following ones:

- Patients under 18 years of age
- Pre-operative complications of the cataract surgery (vitreous humour prolapse, bleeding)
- Expected post-operative problems caused by the cataract surgery
- Previous intraocular operation on the eye to be operated on
- Combined surgical operation at the time of implantation, e.g. glaucoma surgery or keratoplasty; iridectomy is acceptable
- Acute infection or inflammation in the eye

- Chronic use of steroids, immunosuppressant and/or antineoplastic agents
- Significant diseases of the cornea or surface pathology of the eye, e.g. sicca syndrome, keratoconus or blepharitis
- Chronic uveitis, iritis, iridocyclitis or rubeosis iridis
- Iris atrophy
- Aniridia
- Proliferative diabetic retinopathy
- Rubella, congenital/developmental cataracts
- Corneal dystrophy, particularly endothelial dystrophy
- Severe microphthalmia or macrophthalmia
- Decompensated chronic glaucoma
- Systematic diseases with ophthalmological involvement (e.g. diabetes, AIDS complications), that might cause misinterpretation of postoperative findings
- Congenital cataracts in both eyes
- Patients with a history of recurrent anterior or posterior segment inflammations of unknown etiology
- Patients with shallow anterior segment, e.g. microphthalmia or certain types of chronic angle closure glaucoma
- Patients in whom an intraocular lens might interfere with the ability to evaluate, diagnose or treat posterior segment diseases
- Implantation in the functionally better eye, or in patients having only one eye

Certain operative complications may lead to a contraindication against the use of the device, e.g.:

- Rupture of the posterior capsular bag
- Detachment of Descemet's membrane
- Heavy anterior chamber bleeding
- Continued bleeding
- Iris damage
- Vitreous loss
- Irremovable vitreous bulge in the anterior chamber
- Uncontrollable intraocular pressure

COMPLICATIONS OF THE IMPLANTATION (SPECIFIC AND NON-SPECIFIC)

- Shallow anterior chamber
- Corneal oedema
- Corneal dystrophy
- Pupillary block
- Iridocyclitis and vitritis
- Uveitis
- Precipitates on the intraocular lens
- Decentration of the intraocular lens
- (Sub-) luxation of the intraocular lens
- Infection and panophthalmitis

INSTRUCTIONS FOR USE

- Hyphema
- Glaucoma
- Vitreous bulge
- Vitreous detachment
- Retinal detachment
- Post-cataract membrane / Secondary membrane (after cataract)
- Inflammatory pupillary membrane
- Iris atrophy
- Cystoid macular oedema
- Anterior capsular contraction syndrome
- Undesired ametropia
- Aniseikonia

SIDE EFFECTS

There have been occasional reports of secondary glaucoma after the implantation of posterior chambers lenses in patients with controlled glaucoma.

Occasionally, as a result of intra-operative complications of the cataract extraction, post-operative problems such as macular edema, hyphema or vitritis occur, so such patients should be carefully monitored post-operatively.

WARNINGS

- Do not use if packaging is open or damaged.
- Rinse the implant only with sterile rinsing solutions such as sterile Ringer's solution or sterile BSS®.
- The implant cannot be re-sterilised.
- Do not use the implant if the expiration date has expired.
- Do not re-use the implant.
- Re-use can lead to reduced effectiveness or functionality of the implant and/or loss of sight.
- An implant that is damaged during or by an implantation must be disposed immediately.
- Improper handling of the implant may cause injury.
- Improper handling and/or folding may lead to damage of the implant.
- Consult accompanying documents and warning instructions.
- Consult Instructions for Use.
- Do not use implants from dried-out packages.
- The implant must remain moistened while exposed to the air.

PRECAUTIONARY MEASURES

- The implant must be used by an ophthalmologist/ophthalmic surgeon only.
- Before performing an implantation, the surgeon must carefully read the Instructions for Use and all the information material provided by MORCHER® GmbH for the correct handling and insertion of the implant.

STERILISATION AND PACKAGING

- The implant is sterilised by pressurized steam (autoclav).
- The implant is supplied wet and sterile.
- The primary container is sealed, germproof, and packed into a sealed and germproof transparent pouch (sterile double packaging).
- Sterility is only guaranteed as long as the transparent pouch and the primary container are neither opened nor damaged.
- The expiration date is printed on the folded box and the primary container.

STORAGE CONDITIONS / TRANSPORT CONDITIONS

- Store at room temperature +15 °C – +25 °C
- Transport at -18 °C – +60 °C
- Protect from direct sunlight and extreme temperatures
- Protect from moisture / water and store in a dry environment

PATIENT INFORMATION

- Each patient should be fully informed about the implant before implantation.
- The surgeon performing the implantation must inform the patient of all side effects and risks, particularly those that have not already been pointed out in the initial consultation.
- The patient must be instructed to inform the doctor about all side-effects, especially those to which attention has not already been drawn during the initial consultation.

PATIENT CARD

The relevant details and one label must be attached onto the enclosed patient card. The date of the implantation must be written on this card and the eye that was operated on must be crossed off. This card must be given to the patient for safekeeping, and it should be presented to every eye specialist whom the patient consults.

NOTES ON HANDLING / IMPLANTATION

1. Check all labels with regard to type, specific details and expiration date before opening the transparent pouch.
2. Only pick up the implant with smooth-tipped forceps and handle it carefully to avoid damage.
3. Under aseptic conditions, remove the implant as follows:
 - Open the transparent pouch at the points indicated and remove the primary container
 - Carefully pull off the seal from the primary container
 - Carefully remove the implant; Do not pull on the implant.
4. Check implant before use and only use an undamaged implant.
5. Prior to the implantation it is recommended to rinse the implant with sterile Ringer's solution or sterile BSS® - solution.

INSTRUCTIONS FOR USE

Folding the IOL placed in the inlet (Graphics can be found on pages 4-5)

1. The IOL is placed in a sealed liquid-filled container and is kept by a folding aid.
2. The lens is fixated in the folding aid with the curve facing downwards and the haptic angulation facing upwards.
3. Keep the IOL moistened during the folding process.
4. Press the longitudinal sides of the folding aid gently together to fold and remove the lens by using a thin forceps with smooth and round ends. If you want to fold the implant in a different way or not at all, press the folding aid gently together on the short side. The lens will be released and can be removed.
5. The lens will gently unfold when it is rotated in the capsular bag.

Folding of the IOL without inlet

1. It is recommended to fold the lens along the 6-12 o'clock axis. Do not deviate more than 10 degrees from this axis. If the lens is folded along another axis, this could result in damage to the haptics-optics attachment, to tears in the optics or a break in the haptics.
2. It is important to coat the lens with sufficient viscoelastic before folding with forceps.
3. Lens design and power must be considered when selecting the forceps, and the incision size must be considered accordingly.

After implantation of the IOL, please ensure that the IOL is completely unfolded and correctly positioned.

IMPORTANT NOTES

The implant has not been clinically investigated in pregnant or nursing women. Ophthalmologists/ophthalmic surgeons considering lens implantation in such patients should carefully assess the benefit-risk-ratio before implantation.

A positive outcome of intraocular lens implantation cannot be assured in patients with pre-existing ocular disorders (chronic medicinal miosis, diabetic retinopathy, previous cornea transplant, history of retinal detachment or iritis, glaucoma, corneal disease etc.). Ophthalmologists/ophthalmic surgeons considering lens implants in such patients should explore the use of alternative methods. Close postoperative IOP monitoring is indicated in implant patients with pre-existing ocular disorder.

It should be kept in mind that central visual acuity may not be improved in patients with age-related macular or retinal degeneration.

A secondary iridectomy for pupillary block may be avoided if one or more iridectomies are performed at the time of intraocular lens implantation.

An IOL with a small optic may result in a patient experiencing glare under certain light conditions with large pupils. Ophthalmologists/ophthalmic surgeons should consider this before implanting an IOL having a small optic.

Any ingestion and/or application of phosphate-containing substances before/during/after surgery must be avoided, as this can lead to calcification of the implant and result in loss of vision.

In the case of patients with an impaired calcium-phosphate balance (e.g. patients with diabetes), there is an increased risk of the IOL possibly becoming more cloudy at a later point (calcification).

Colouring substances in connection with hydrophilic implants should be avoided whereas this may cause permanent colouring of the implant material.

Notes for laser treatment:

Focus the laser beam precisely on the action site behind the implant. A laser beam focused on the implant itself will lead to a damage of the implant.

DISPOSAL

- The packaging and its contents must be disposed in accordance with the current local regulations.
- There are no special disposal regulations applicable for products which are uncontaminated.
- Contaminated products / packaging must be handled and disposed in accordance with the provisions of the respective clinic.

EXCLUSION OF LIABILITY

MORCHER® GmbH cannot be held liable for the surgical techniques used by the surgeon, nor for the selection of the implant for any eye. This is solely the responsibility of the surgeon.

RETURN OF A DAMAGED/NON-STERILE IMPLANT

If an implant is returned, the serial number, batch code and reason for its return must be provided. The packaging containing the returned item must be in accordance to the condition of the implant and, where necessary, must warn of any contamination.

INSTRUCTIONS FOR USE

REPORTING ON ADVERSE REACTIONS

Adverse reactions and/or complications affecting vision, which can be attributed to the implant and were not previously expected with regard to nature, severity or degree of the incidence, must be reported to MORCHER® GmbH. This information is requested from all surgeons, in order to document long-term effects caused by use of the implant in a timely manner and, if possible, to eliminate them. In such cases, please contact:

MORCHER® GmbH

Quality Assurance Department
Kapuzinerweg 12
70374 Stuttgart
Germany

Phone	+49 (0) 711 / 95 320-0
Fax	+49 (0) 711 / 95 320-80
E-Mail	complaints@morcher.com

ABBREVIATIONS AND PACKAGING SYMBOLS

ABBREVIATION / SYMBOL	MEANING
IOL	Intraocular lens
CTR	Capsular Tension Ring
PC	Posterior Chamber
AC	Anterior Chamber
ACCC	Anterior Continuous Curvilinear Capsulorhexis
PCCC	Posterior Continuous Curvilinear Capsulorhexis
ICCE	Intracapsular Cataract Extraction
ECCE	Extracapsular Cataract Extraction
PCO	Posterior Chamber Opacification
D	Diopter in aqueous humour
ØB	Body diameter
ØT	Total diameter
ØO	Optic diameter
ØI	Inner diameter
R	Radius
Th	Haptic side view
A - CON	A - constant
AC DEPTH / ACD	Anterior Chamber Depth
PMMA	Polymethylmethacrylate
PEMA	Polyethylmethacrylate
UV	Ultraviolet

INSTRUCTIONS FOR USE

HEV	High Energy Visible light
C-F-M	Compression forged method
CQ	Clinical Quality
BSS®	Balance Salt Solution
SN	Serial Number
UDI	Unique Device Identification
GTIN	Global Trade Item Number
MD	Medical Device
REF	Catalogue number
LOT	Batch code
	Use by: Year - Month
 EXP	Use by: Year - Month
	Manufacturer
	Country and date of manufacture
	Importer
	Distributor
CE XXXX	CE Mark with 4 - digit number of the notified body
STERILE R	Sterilised using irradiations
STERILE 	Sterilised using steam
STERILE P	Sterilised using plasma
	Double sterile barrier system
kGy	Kilogray
	Do not re - use

Do not re-sterilise

Consult Instructions for Use

Caution, consult accompanying documents

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