INSTRUCTIONS FOR USE

EYEJET® CTR

SINGLE USE PRE-LOADED CAPSULAR TENSION RING

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


INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

The MORCHER® EyeJet® CTR consists of a MORCHER® Capsular Ring (CR) and a single-use injector MORCHER® EyeJet®, which already contains the preloaded Capsular Ring.

MORCHER® Capsular Ring, hereinafter also referred to as implant, is a sterile, non-optical medical device designed to stabilize the capsular bag after cataract extraction. MORCHER® Capsular Rings are made of polymethylmethacrylate (PMMA), and are either blue in colour or colourless.

The MORCHER® single-use injector MORCHER® EyeJet® CTR, hereinafter referred as injector, is an accessory that facilitates the implantation of the MORCHER® Capsular Ring.

Once the preloaded MORCHER® EyeJet® CTR is unpacked, the implant is already hooked in the injector and should be immediately drew into the injector completely. Subsequently, the injector is ready for use in a surgery for implanting of a CR and an intraocular lens (IOL).

There are two groups of pre-loaded MORCHER® Capsular Rings:

- Capsular Tension Ring (CTR)
- Cionni Ring

The MORCHER® EyeJet® CTR is available in two different models:

1. For clockwise implantation (RIGHT)
2. For anti-clockwise implantation (LEFT)

Graphic can be found on pages 4 + 5.

INTENDED USE

MORCHER® EyeJet® CTR is intended for the stabilisation of the crystalline lens capsule in the presence of weak or partially absent zonules in patients undergoing cataract extraction with intraocular lens implantation.

INDICATIONS

- Used only in connection with a cataract operation with IOL.
- Only in connection with IOL's conform to the international standard Ophthalmic implants – Intraocular lenses – (ISO 11979 series)
- MORCHER® Capsular Rings help to stabilize the capsular bag in the case of high myopia.

Additional indications only for Capsular Tension Rings and Cionni Rings:

- With Zonulolysis
- In Pseudoexfoliation
- Primary zonular weakness
- Defective or missing zonules

Additional indications only for Capsular Tension Rings:

- Marfan's syndrome

Additional indications only for Cionni Rings:

- Secondary zonular weakness (e.g. Marfan's syndrome)

CONTRAINDICATIONS

The usual contraindications for a cataract operation apply here, as well as the following ones:

- Patients, whose eyes are not fully developed, i.e. implants are contraindicated in eyes that are still in their developmental stage.
- Pre-operative complications of the cataract surgery (vitreous humour prolapse, bleeding)
- Expected post-operative problems caused by the cataract surgery
- Severe microphthalmia or macrophthalmia
- Patients with a flat front section, e.g. microphthalmia or certain types of chronic angle closure glaucoma
- Chronic uveitis
- ICCE (intracapsular cataract extraction) technology
- Flat anterior chamber
- Progressive ophthalmic diseases (proliferative diabetic retinopathy, uncontrolled glaucoma)

Certain operative complications could also lead to a contraindication against the use of the implant, e.g.:

- Rupture of the capsule of the lens with or without vitreous prolapsed.
- Damage to the capsular bag
- Prolonged bleeding or other visual impairments
- Anterior chamber inflammation (occurring in nearly half of patients)
- Cortical remnants
- Vitreous problems
- Druses
- Cortical macular edema
- Macular degeneration
- Posterior capsular opacification (at least in conventional capsular rings without sharp edge)
- IOL decentered
- Elevated IOP
- Optic atrophy
- Retinal detachment
- Deposits on IOL
- Iritis
- Cystoid macular edema
- Posterior synechiae

COMPLICATIONS OF THE IMPLANTATION (SPECIFIC AND NON-SPECIFIC)

- angle closure

SIDE EFFECTS

The possible side effects of the Capsular Rings are as follows:

- Zonular trauma
- Impairments caused during cortical cleansing
- Possible dislocation of the IOL or the capsular bag when the posterior capsular surface tears

WARNINGS

- The implant and the EyeJet® CTR may not be used if the sterile packaging is already open or damaged.
- The implant and the injector cannot be re-sterilized.
- The implant and the injector are only suitable for single use.
- The implant may only be rinsed with sterile intra-ocular rinse solutions such as sterile Ringer solution or sterile BSS®-solution.
- The injector may only be used for the implant supplied with it.
- The implant may not be re-used. Re-use can lead to reduced effectiveness or the loss of the implant's functionality. Moreover, it can cause loss of sight in the eye.
- The injector may not be re-used because it is only designed for single-use and its materials are not suitable for reuse/reprocessing.
- The accompanying documents and the warning instructions must be observed.
- The directions for use must be consulted.
- Eyes with pseudoexfoliation syndrome and flat anterior chamber run a major risk of zonular defects or other intraoperative complications.
- An implant that has been damaged during implantation should not be allowed to remain in the eye!
- If a posterior capsulorhexis (PCCC) is created, the MORCHER® Capsular Ring must be used after the PCCC. Otherwise there is a risk that the tension of the MORCHER® Capsular Ring may cause the capsulorhexis to lacerate uncontrollably.
- Warning for Cionni Rings: In the case of Cionni Rings for which MORCHER® recommends an injector, it is important to note that when retracting you stop at the right time so that the fixation loops do not get damaged by striking against the cannula. The fixation loops remain outside the cannula.
- Do not use the implant if the expiration date has passed.

INSTRUCTIONS FOR USE

PRECAUTIONARY MEASURES

- The surgeon performing the implantation must carefully read all the information material provided by MORCHER® GmbH for the proper use and insertion of the implant before carrying out the implantation.
- The implant may only be used by a doctor experienced in the implantation of intraocular lenses, or under his/her supervision. The use of this product requires a detailed understanding of the technical fundamental principles, clinical applications and the risks of this product.

Note when using lasers

The laser beam must be focused directly on the target point behind the implant. If a laser beam is focused on the implant, it can result in damage to the latter.

STERILISATION AND PACKAGING

- The MORCHER® EyeJet® CTR is sterilised by gamma irradiation.
- The MORCHER® EyeJet® CTR is supplied dry and sterile.
- The primary blister is sealed, germproof, and packed into a sealed and germproof secondary blister (sterile double packaging).
- Sterility is only guaranteed as long as both blisters are neither opened nor damaged.
- The expiration date is printed on the folded box and the primary container.

STORAGE CONDITIONS/TRANSPORT CONDITIONS

- Store at room temperature +15 °C – +25 °C
- Transport at -18 °C – +40 °C
- Protect from direct sunlight and extreme temperatures
- Protect from moisture / water and store in a dry environment

PATIENT INFORMATION

- Each patient should be fully informed about the implant before implantation.
- The surgeon performing the implantation must inform the patient of all side effects and risks, particularly those that have not already been pointed out in the initial consultation.
- The patient must be instructed to inform the doctor about all side-effects, especially those to which attention has not already been drawn during the initial consultation.

PATIENT CARD

The relevant details and one label must be attached onto the enclosed patient card. The date of the implantation must be written on this card and the eye that was operated on must be crossed off. This card must be given to the patient for safekeeping, and it should be presented to every eye specialist whom the patient consults.

NOTES ON HANDLING/IMPLANTATION

1. Check all labels with regard to type, specific details and expiration date before opening the blister packs of the MORCHER® EyeJet® CTR.
2. Under aseptic conditions, remove the device as follows:
 - Open the outer blister pack at the points indicated and remove the primary container
 - Carefully pull back the seal of the second blister pack
 - Carefully remove the MORCHER® EyeJet® CTR containing the preloaded CR.
3. Check the device before use and only use an undamaged device.
4. The correct handling of the MORCHER® EyeJet® CTR is shown in diagrams on pages 4 + 5.
 - The preloaded CR should be drawn into the injector completely and the clamp of the injector should be removed. The injector is now ready for use.
5. Implantation:
 - Before implantation, fill the capsular bag with sufficient viscoelastic.
 - During the implantation, please ensure that the writing on the injector is facing upwards. This guarantees that the injector has the right position to eject the implant. Slowly push the tip of the injector through the capsulorhexis until it reaches approx. the middle of the capsular bag.
 - Push the plunger carefully to implant the MORCHER® CR in the capsular bag.
 - The implantation must be performed using a direction of feed that is tangential to the extent possible to prevent any probable interlocking of the leading eyelets. The implants are designed in such a manner that it follow the natural curve of the capsular bag.
 - If required, the implant must be drawn backwards and the fornix of the capsular bag distended by injecting a highly viscous viscoelastic agents.

- Ensure that the plunger has been fully pushed through to its end position.
- Carefully remove the hook of the injector from the eyelet of the CR and pull the hook back into the injector so that no patient's tissue can be damaged when removing the injector.
- If a posterior capsulorhexis (PCCC) is intended by the surgeon, the MORCHER® CR must be used after the PCCC. Otherwise there is a risk that the tension of the MORCHER® Capsular Ring may cause the capsulorhexis to lacerate uncontrollably.

Instruction for immediate implantation of an intraocular lens

The cradle of the implanted open loop lens or the axis of the intraocular lens must be set at a 90° angle to the axis of the capsular ring. It is an advantage to have the cradle of the lens positioned posterior to the ring in the fornix of the capsular bag. This can be achieved by applying pressure on the lens in the region of the haptic or on the cradle itself.

Additional instructions for Cionni Ring

The size of the anterior rhexis depends on the implant used. Corresponding information can be requested from MORCHER® GmbH.

When inserting a Cionni Ring, ensure that the leading fixation loop (the curved end of the ring) with the eyelet is located outside of the capsular bag.

Position the fixation eyelets of the Cionni Ring at the position of the zonular defect.

Once the Cionni Ring is in the correct position, it should be sutured to the sclera using the fixation eyelet.

A double-armed 9-0 Prolene suture should be used for fixation of the implant.

IMPORTANT NOTES

The implants have not been clinically investigated in pregnant or nursing women. Ophthalmologists/ophthalmic surgeons considering lens implantation in such patients should carefully assess the benefit-risk-ratio before implantation.

DISPOSAL

- The packaging and its contents must be disposed in accordance with the current local regulations.
- There are no special disposal regulations applicable for products which are uncontaminated.
- Contaminated products/packaging must be handled and disposed in accordance with the provisions of the respective clinic.

EXCLUSION OF LIABILITY

MORCHER® GmbH cannot be held liable for the surgical techniques used by the surgeon, nor for the selection of the implant for any eye. This is solely the responsibility of the surgeon.

RETURN OF A DAMAGED/NON-STERILE IMPLANT

If an implant is returned, the serial number, batch code and reason for its return must be provided. The packaging containing the returned item must be in accordance to the condition of the implant and, where necessary, must warn of any contamination.

REPORTING ON ADVERSE REACTIONS

Adverse reactions and/or complications affecting vision, which can be attributed to the implant and were not previously expected with regard to nature, severity or degree of the incidence, must be reported to MORCHER® GmbH. This information is requested from all surgeons, in order to document long-term effects caused by use of the implant in a timely manner and, if possible, to eliminate them. In such cases, please contact:

MORCHER® GmbH

Quality Assurance Department
Kapuzinerweg 12
70374 Stuttgart
Germany

Phone	+49 (0) 711/95 320-0
Fax	+49 (0) 711/95 320-80
E-Mail	complaints@morchher.com

INSTRUCTIONS FOR USE

ABBREVIATIONS AND PACKAGING SYMBOLS

ABBREVIATION / SYMBOL	MEANING
IOL	Intraocular lens
CTR	Capsular Tension Ring
PC	Posterior Chamber
PCCC	Posterior Continuous Curvilinear Capsulorhexis
ICCE	Intracapsular Cataract Extraction
ØB	Body diameter
ØT	Total diameter
ØO	Optic diameter
ØI	Inner diameter
R	Radius
Th	Haptic side view
PMMA	Polymethylmethacrylate
C - F - M	Compression forged method
CQ	Clinical Quality
BSS®	Balanced Salt Solution
SN	Serial Number
	Unique Device Identification
GTIN	Global Trade Item Number
	Medical Device
	Catalogue number
	Batch code
	Use by: Year - Month
 EXP	Use by: Year - Month
	Manufacturer
	Country and date of manufacture
	Importer
	Distributor
	CE Mark with 4-digit number of the notified body
	Sterilised using irradiations

	Double sterile barrier system
kGy	Kilogray
	Do not re-use
	Do not re-sterilise
	Consult Instructions for Use
	Caution, consult accompanying documents
	Keep dry
	Keep away from sunlight
	Temperature limit
	Do not use if packaging is damaged
	This side up
	Recyclable packaging

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