Insturctions for Use  
Vitreous Cutter

The MID Labs Vitreous Cutter is intended for use/intended purpose to remove vitreous and intra-ocular tissue, when used in conjunction with any compatible ophthalmic surgical system.

System Requirements:				STERILE R	
REF	Minimum Actuation Pressure	Connection Type	Probe Type		
2350CE / 2370CE / 2390CE	138 KPa (20.0 psi)	Luer (male)	Conventional		
2370BB / 2390BB / 2310BB	138 KPa (20.0 psi)	Luer (male)	Bi-Blade		
2450CE / 2470CE / 2490CE	172 KPa (25.0 psi)	Luer (male)	Conventional		
2470BB / 2490BB / 2410BB	172 KPa (25.0 psi)	Luer (male)	Bi-Blade		
2470QB / 2490QB / 2410QB	172 KPa (25.0 psi)	Luer (male)	Bi-Blade		
2570BB / 2590BB / 2510BB	200 KPa (29.0 psi)	Luer (male)	Bi-Blade		
2670BB / 2690BB / 2610BB	214 KPa (31.0 psi)	CPC (male)	Bi-Blade		
2870BB / 2890BB / 2810BB	259 KPa (37.5 psi)	CPC (male)	Bi-Blade		
2970BB / 2990BB / 2910BB	190 KPa (27.5 psi)	CPC (female)	Bi-Blade		

- INSTALLATION**
1. Tubing set is color coded for easy identification; GREEN striped line for cutting actuation pressure and CLEAR tubing for aspiration.
  2. Handle using appropriate aseptic technique.
  3. Attach GREEN striped tubing to the vitrectomy actuation connector on system console.
  4. Attach CLEAR tubing to the aspiration connector of the phaco tubing set.

**TESTING AND PRIMING**

*Note: Refer to user manual of system console for additional instruction on priming and testing the cutter.*

1. Place needle of cutter in a container of sterile solution and actuate cutting. Look for bubbles or other leaks. If bubbles are seen at cutting port, do not use cutter.
2. Under the microscope, test for appropriate cutting action of the vitreous cutter. Do not run the cutter outside sterile solution for more than a few moments. If the vitreous cutter is not operating properly, do not use cutter.
3. Place needle of cutter in a container of sterile solution and prime tubing set. To prime, actuate aspiration until all bubbles are removed from tubing set.

Glossary of Symbols							
Symbol	Title/Description	Applicable Standard	Reference Number	Symbol	Title/Description	Applicable Standard	Reference Number
	Manufacturer	EN ISO 16223-1	6.1.1		Do not re-use	EN ISO 16223-1	6.4.2
	Authorized Representative in the European Community	EN ISO 16223-1	6.1.2		Do not use if package is damaged	EN ISO 16223-1	6.2.8
	Use-by date	EN ISO 16223-1	6.1.4		Refer to instruction manual/booklet	ISO 7010	M002
	Batch code	EN ISO 16223-1	6.1.6		Caution	EN ISO 16223-1	6.4.4
	Catalogue number	EN ISO 16223-1	6.1.6		Caution: Federal (U.S.) Law Restricts this device to sale by or on the order of a physician.	21 CFR 801	§ 801.109(b)(1)
	Sterilized using irradiation	EN ISO 16223-1	6.2.4		European Conformity - CE Mark	Council Directive 93/42/EEC	ANNEX XII
	Do not re-sterilize	EN ISO 16223-1	6.2.6		Waste electrical and electronic equipment: Refer to local regulations for disposal	WEEE Directive 2012/19/EU	ANNEX IX
	Medical Device	EN ISO 16223-1	6.7.7		Caution, sharp edges	IEC 60417	6043
	Double sterile barrier system	EN ISO 16223-1	6.2.12		Translation	EN ISO 16223-1	6.7.8
	Vitrectomy	BS EN 80801-2-48	ANNEX D		Contains or present of Phthalates (DEHP)	EN 16986	4.2
	Authorized Representative in Switzerland	Medical Devices Ordinance (MedDO)	Art. 4 para. 1 let. g				

If any MID Labs products are received in a defective condition, your distributor or MID Labs is to be notified immediately.

**Medical Instrument Development Laboratories, Inc.**  
557 McCormick Street, San Leandro, CA 94577 U.S.A.  
Phone: 1-800-929-5227 (inside U.S. only) Phone: 510-357-3952  
Fax: 510-357-1582 Email: info@midlabs.com

Each MID Labs product is identified by a lot number. Refer to this number when discussing the product with your MID Labs representative.

Use of accessories other than those of MID Labs may affect system performance and create potential hazards.

MID Labs assumes no responsibility for complications that may arise from use of accessories other than those of MID Labs or improper cleaning, sterilization, or use, of MID Labs products.

Any serious incident that has occurred in relation to this device should be reported to MID Labs and the competent authority of the Member State in which the user and/or patient is established.

EC

REP

EMERGO EUROPE  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands

0197

CH

REP

MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland

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MI Medical Instrument Development Laboratories, Inc