Instructions for Use

BioComFold

Foldable Intraocular Lenses

Morchers GmbH

Kapuzinerweg 12 | 70374 Stuttgart | Germany

www.morcher.com

CE 0483

..... 3

[en] ENGLISH | INSTRUCTIONS FOR USE

BioComFold | Foldable Intraocular Lenses 5

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





DISCLAIMER

Please note: Not all models are available in all countries and from all distributors.

PRODUCT DESCRIPTION

The BioComFold (BCF) intraocular lenses (IOL), hereafter referred to as lenses or implants, are sterile optical medical implants designed to replace the natural lens. The quality of the material used permits folding the optic for small incision insertion. The lenses are made from one piece of a hydrophilic acrylic copolymer including a copolymerized UV filter and, if applicable, an additional copolymerized HEV filter (high energy visible light filter / blue light blocker). Each lens is tested separately for its optical quality and refractive power.

BCF IOL are available with different diameters, different sizes of optic and in various forms in order to ensure optimal fitting of the implant in the eye.

DEVICE MATERIALS

The BCF IOL are made from a hydrophilic acrylic copolymer. This polymer contains < 1 % of a copolymerizable UV absorber and, if applicable, additionally < 0.3 % of a copolymerizable HEV filter (blue light blocker).

SPECIFICATIONS

Model (REF)	Total Ø [mm]	Body Ø [mm]	Remarks
CARE FLEX 2	11.0	6.0	---
CARE FLEX2 A	11.0	6.0	aspheric design
41B	10.0	5.5	---
46G	15.0	6.0	includes HEV filter
90L	15.0	6.5	---
90S	15.0	6.4	---
92B	15.0	7.0	---
92C	10.5	6.0	---
92G	15.0	6.5	---
92S	15.0	6.0	---
95S	15.0	7.0	---
97B	10.5 – 11.0	6.0 – 6.5	---
97L	10.5	6.0	---

Refractive index = 1.46 for all BCF IOL.

INTENDED PURPOSE

The BioComFold (BCF) intraocular lenses are intended to serve as replacement of the natural lens. The lenses are available for fixation into the capsular bag. They are used after ECCE (extra-capsular cataract extraction).

INTENDED USER

BCF IOL are used by either an ophthalmic surgeon or an ophthalmologist during a state-of-the-art operation.

USE ENVIRONMENT

BCF IOL are used under aseptic conditions in an operating room or a doctor's practice using a surgical microscope.

TARGET GROUP

Target group are adults 18 years and older.

The implants have not been clinically investigated in pregnant or nursing women. Ophthalmologists/ ophthalmic surgeons, considering implantation in such patients, should carefully assess the benefit-risk ratio before implantation.

INDICATIONS

- For senile cataracts
- As correction of an aphakic eye where other optical aids are not available, or are not implementable
- For traumatic cataracts
- In exceptional situations with congenital or early acquired cataracts

CONTRAINDICATIONS

The usual contraindications for a cataract operation apply here, as well as the following ones:

- Patients under 18 years of age
- Pre-operative complications of the cataract surgery (vitreous prolapse, hemorrhage)
- Anticipated post-operative problems of cataract surgery
- Previous intraocular operation on the eye scheduled for surgery
- Combined surgical intervention at the time of implantation, e.g. glaucoma surgery or keratoplasty (iridectomy is acceptable)
- Acute infection or inflammation of the eye
- Chronic use of steroids, immunosuppressant and/or antineoplastic agents
- Significant corneal diseases or surface pathology of the eye, e.g. sicca syndrome, keratoconus or blepharitis
- Chronic uveitis, iritis, iridocyclitis or rubeosis iridis
- Iris atrophy
- Aniridia

- Proliferative diabetic retinopathy
- Rubella, congenital/developmental cataracts
- Corneal dystrophy, particularly endothelial dystrophy
- Severe microphthalmia or macrophthalmia
- Decompensated chronic glaucoma
- Systemic diseases with ophthalmological involvement (e.g. diabetes, complications of AIDS) that might cause misinterpretation of postoperative findings
- Congenital cataracts in both eyes
- Patients with a history of recurrent anterior or posterior segment inflammations of unknown etiology
- Patients with a shallow anterior segment, e.g. microphthalmia or certain types of chronic angle closure glaucoma
- Patients in whom an intraocular lens might interfere with the ability to evaluate, diagnose or treat posterior segment diseases
- Implantation in the functionally better eye, or in patients having only one eye

Certain operative complications may lead to a contraindication against the use of the device, such as:

- Rupture of the posterior capsular bag
- Detachment of Descemet's membrane
- Significant anterior chamber bleeding
- Persistent bleeding
- Iris damage
- Vitreous loss
- Irremovable vitreous bulge in the anterior chamber
- Uncontrollable intraocular pressure

INTENDED CLINICAL BENEFIT

The BCF IOL are used to replace the natural lens in case of cataract or aphakia. They lead to the benefit of a restored or improved visual acuity. Due to the material properties, the lenses can be folded and implanted through a small incision.

PERFORMANCE CHARACTERISTICS

BCF IOL significantly improve visual acuity as demonstrated by scientific literature and post-market clinical follow-up (PMCF) activities of the manufacturer.

It fulfils the requirements of the DIN EN ISO 11979-7 on clinical investigations on IOL as proven by a PMCF (post-market clinical follow-up) study on an open-loop BCF IOL: a postoperative decimal visual acuity (VA) of 0.5 or better was achieved in 96.3 % (52/54) of cases without ocular comorbidity and in 91.0 % (131/144) of all cases at last regular visit (3 to 14 months postoperatively).

Data reported in literature on BCF IOL and similar devices as well as data collected by PMCF questionnaires confirm these data.

PMCF questionnaires collected data on surgery/short-term follow-up on > 10,000 cases, data on long-term follow-up (3 months to > 10 years) on > 4,900 cases. Post-operative improvement of VA was reported for 70 % to 100 % of cases; a decimal Best-corrected VA of ≥ 0.5 (≤ 0.3 logMAR) was achieved in 75 % to 100 % of the cases. VA was considered clinically appropriate by the ophthalmologist/ophthalmic surgeon in 90 - 100 % of the cases.

For details see section “Summary of safety and clinical performance”.

LIFETIME

The lifetime of the BCF IOL is 10 years. Without any complications or adverse events, the device may stay in the patient for his lifetime.

COMPLICATIONS OF THE IMPLANTATION (SPECIFIC AND NON-SPECIFIC)

- Complications associated with the general surgical procedure
- Shallow anterior chamber
- Corneal oedema
- Corneal dystrophy
- Pupillary block
- Intrapupillary membrane
- Iridocyclitis, vitritis, uveitis
- Precipitates on the intraocular lens (calcification/opacification of the intraocular lens)
- Decentration of the intraocular lens
- (Sub-) luxation of the intraocular lens
- Dislocalization of the intraocular lens
- Iris atrophy
- Infection, endophthalmitis and panophthalmitis
- Suprachoroidal hemorrhage
- Hyphema
- Anterior synechia
- Uncontrolled elevated intraocular pressure
- Glaucoma
- Vitreous bulge
- Vitreous detachment
- Retinal detachment
- Posterior capsule opacification (PCO)/post-cataract membrane/secondary membrane
- Inflammatory pupillary membrane
- Cystoid macular oedema
- Anterior capsular contraction syndrome (ACCS)
- Capsular block syndrome
- Toxic anterior segment syndrome (TASS)
- Hypopyon

- Undesired ametropia
- Aniseikonia
- Surgically induced astigmatism
- Iris capture
- Dysphotopsia; glare

SIDE EFFECTS

There have been occasional reports of secondary glaucoma after the implantation of posterior chambers lenses in patients with controlled glaucoma.

Occasionally, as a result of intraoperative complications of the cataract extraction, post-operative problems such as macular edema, hyphema or vitritis occur, so such patients should be carefully monitored post-operatively.

WARNINGS

- Do not use if the sterile packaging is already open or damaged.
- Do not re-sterilise the implant. Re-sterilization can lead to damage to the material and loss of functionality of the implant.
- Do not re-use the implant. The implant is intended for single use. Re-use can lead to reduced effectiveness or functionality of the implant and/or loss of sight.
- Do not use the implant if the expiration date has expired.
- Rinse the implant only with sterile rinsing solutions such as sterile Ringer's solution or sterile BSS® (balanced salt solution).
- Do not use a damaged implant. An implant that is damaged during the implantation must be disposed of immediately.
- Improper handling and/or folding of the implant may lead to damage of the implant and may cause injury.
- Observe accompanying documents and warning instructions before using the device.
- Consult Instructions for Use before using the device.
- Do not use implants from dried-out packages.
- The implant must remain moistened while exposed to the air.
- Do not use the implant if it may have been contaminated with unsterile material (e.g. if it falls on the floor). Dispose of the potentially contaminated implant and use a new product.

PRECAUTIONARY MEASURES

QUALIFICATION AND TRAINING FOR USERS

- The implant may only be used by an ophthalmologist/ophthalmic surgeon.
- The implant may only be used by a doctor experienced in the implantation of IOL, or under his/her supervision. The use of the lenses requires a detailed understanding of the technical fundamental principles, clinical applications and the risks of the devices.
- The surgeon performing the implantation should carefully read all the information material provided by the manufacturer for the proper use and insertion of the implant before carrying out the implantation.

CALCIFICATION OF THE IMPLANT

In rare cases opacification/calcification of the IOL may occur, which usually results in loss of vision. The exact pathomechanism has not yet been fully understood. To minimize the risk of IOL calcification, the IOL should not come into contact with phosphate-containing substances. We thus recommend avoiding phosphate-containing substances during surgery.

NOTE WHEN USING LASERS

The laser beam must be focused directly on the target point behind the implant. If a laser beam is focused on the implant, it can result in damage to the implant.

LIMITATIONS

A positive outcome of IOL implantation cannot be assured in patients with pre-existing ocular disorders (chronic medicinal miosis, diabetic retinopathy, previous cornea transplant, history of retinal detachment or iritis, glaucoma, corneal disease etc.). Ophthalmologists/ophthalmic surgeons, considering lens implants in such patients, should explore the use of alternative methods. Close postoperative monitoring of intraocular pressure is indicated in patients with pre-existing ocular disorder.

It should be kept in mind that central visual acuity may not be improved in patients with age-related macular or retinal degeneration.

An IOL with a small optic may result in a patient experiencing glare under certain light conditions with large pupils. Ophthalmologists/ophthalmic surgeons should consider this before implanting an IOL having a small optic.

In the case of patients with an impaired calcium-phosphate balance (e.g. patients with diabetes), there is an increased risk of the IOL possibly becoming more cloudy at a later point (calcification).

Colouring substances in connection with hydrophilic implants should be avoided whereas this may cause permanent colouring of the implant material.

STERILISATION AND PACKAGING

- The implant is sterilised by pressurized steam (autoclave).
- The implant is supplied wet and sterile.
- The primary container is sealed, germproof, and packed into a sealed and germproof transparent pouch (double sterile barrier system).
- Sterility is only guaranteed as long as the transparent pouch and the primary container are neither opened nor damaged.
- The expiration date is printed on the folded box and the primary container.

STORAGE CONDITIONS / TRANSPORT CONDITIONS

- Store flat at room temperature + 15 °C up to + 25 °C.
- Transport at - 18 °C up to + 60 °C.
- Protect from direct sunlight and extreme temperatures.
- Protect from moisture/water and store in a dry environment.

PATIENT INFORMATION

- Each patient shall be fully informed about the implant, especially about potential risks (including any undesirable side effects, warnings, precautions, contraindications, and limitations) and benefits, before implantation.
- The surgeon performing the implantation shall inform the patient of all side effects and risks, particularly those that have not already been pointed out in the initial consultation.
- The patient shall be instructed to inform the doctor about all side effects, which he/she may observe after the implantation, especially in case of malfunction of the device or changes in its performance.
- The patient shall be informed that there are no specific environmental conditions known that can affect the BCF IOL. Medical procedures or examinations like laser treatment, magnetic resonance imaging (MRI) or X ray scan do not negatively impact the performance of the device and do not pose any risk to the device.
- The patient shall be informed that there are no known circumstances, in which the device interfere with certain diagnostic tests, evaluations, or therapeutic treatments or other procedures, or interfere with the equipment used in these procedures. Please note, that implantation is contraindicated in patients, in whom the IOL might interfere with the ability to evaluate, diagnose or treat posterior segment diseases.
- Written information about the devices intended for the patient is available at <https://patient-information.morcher.com>.

IMPLANT CARD

The relevant details and one label shall be attached onto the enclosed implant card. The date of the implantation shall be written on this card and the eye that was operated on shall be crossed off. Instructions for filling out the implant card are available at <https://ifu.morcher.com>.

The completed card shall be given to the patient for safekeeping to ensure traceability of the implant; it should be presented to every eye specialist whom the patient consults.

NOTES ON HANDLING / IMPLANTATION

1. Check all labels with regard to type, specific details and expiration date before opening the transparent pouch.
2. Only pick up the implant with smooth-tipped forceps and handle it carefully to avoid damage.
3. Under aseptic conditions, remove the implant as follows:
 - Open the transparent pouch at the points indicated and remove the primary container.
 - Carefully pull off the seal from the primary container.
 - Carefully remove the implant. Do not pull on the implant.
4. Check implant before use and only use an undamaged implant.
5. Prior to the implantation it is recommended to rinse the implant with sterile Ringer's solution or sterile BSS®.

FOLDING THE IOL PLACED IN THE INLET (GRAPHICS CAN BE FOUND ON PAGES 3 + 4)

1. The IOL is placed in a sealed liquid-filled container and is kept by a folding aid.
2. The lens is fixated in the folding aid with the curve facing downwards and the haptic angulation facing upwards.
3. Keep the IOL moistened during the folding process.
4. Press the longitudinal sides of the folding aid gently together to fold and remove the lens by using a thin forceps with smooth and round ends.

If you want to fold the implant in a different way or not at all, press the folding aid gently together on the short side. The lens will be released and can be removed.

5. The lens will gently unfold when it is rotated in the capsular bag.

FOLDING OF THE IOL WITHOUT INLET

1. It is recommended to fold the lens along the 6 - 12 o'clock axis. Do not deviate more than 10 degrees from this axis. If the lens is folded along another axis, this could result in damage to the haptics-optics attachment, to tears in the optics or a break in the haptics.
2. It is important to coat the lens with sufficient viscoelastic before folding with forceps.
3. Lens design and power must be considered when selecting the forceps, and the incision size must be considered accordingly.

After implantation of the IOL, please ensure that the IOL is completely unfolded and correctly positioned.

A secondary iridectomy for pupillary block may be avoided if one or more iridectomies are performed at the time of IOL implantation.

USE WITH OTHER DEVICES AND / OR GENERAL PURPOSE EQUIPMENT

The lenses can be used with the accessory MORCHER® Caliper Ring. The Caliper Ring is an accessory designed to size and centre the position of the anterior capsulorhexis. Caliper Rings are available in different sizes.

There are no restrictions known regarding the combination of the lenses with general surgical equipment, ophthalmic viscosurgical devices (OVD, viscoelastics), or intraocular physiological rinsing solutions.

To minimize the risk of IOL calcification, phosphate-containing substances should be avoided during surgery.

DISPOSAL

- The packaging and its contents shall be disposed of in accordance with the current local regulations.
- There are no special disposal regulations applicable for devices, which are uncontaminated.
- Contaminated devices/packaging shall be handled and disposed of in accordance with the provisions of the respective clinic.

EXCLUSION OF LIABILITY

Morcher GmbH and its distributors cannot be held liable for any procedures, methods, or surgical techniques used by the surgeon, nor for the selection of the implant for any eye. This is solely the responsibility of the surgeon.

COMPLAINTS

Complaints should be addressed to

Morcher GmbH

Quality Assurance Department
Kapuzinerweg 12
70374 Stuttgart
Germany

Phone + 49 (0) 711/95 320 - 0
Fax + 49 (0) 711/95 320 - 80
E-Mail complaints@morcher.com

or to your distributor.

REPORTING ON ADVERSE REACTIONS

Any serious incident that occurred in relation to the device should be reported to

Morcher GmbH

Quality Assurance Department
Kapuzinerweg 12
70374 Stuttgart
Germany

Phone + 49 (0) 711/95 320 - 0
Fax + 49 (0) 711/95 320 - 80
E-Mail complaints@morcher.com

and to the competent authority of the member state in which the user and/or patient is established.
Please also inform your distributor if necessary.

RETURN OF A DAMAGED / NON-STERILE DEVICE

If an implant is returned, the serial number, batch code and reason for its return must be provided.
The packaging, containing the returned item, must be in accordance to the condition of the implant
and, where necessary, must warn of any contamination.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The summary of safety and clinical performance (SSCP) referred to in Article 32 of the Medical Device Regulation (EU) 2017/745 is linked to the Basic UDI - DI number 4048509IOL01S3 and is available at <https://ifu.morcher.com> or in the European database on medical devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed> (as soon as this EUDAMED functionality is available).

ABBREVIATIONS AND PACKAGING SYMBOLS

ABBREVIATION / SYMBOL	DESCRIPTION
IOL	Intraocular Lens
CTR	Capsular Tension Ring
PC	Posterior Chamber
AC	Anterior Chamber
ACCC	Anterior Continuous Curvilinear Capsulorhexis
PCCC	Posterior Continuous Curvilinear Capsulorhexis
ICCE	Intracapsular Cataract Extraction
ECCE	Extracapsular Cataract Extraction
PCO	Posterior Chamber Opacification
D	Diopter in aqueous humour
ØB	Body diameter
ØT	Total diameter
ØO	Optic diameter
ØI	Inner diameter
R	Radius
Th	Haptic side view
A-CON	A-constant
AC DEPTH/ACD	Anterior Chamber Depth
PMMA	Polymethylmethacrylate
PEMA	Polyethylmethacrylate
UV	Ultraviolet
HEV	High Energy Visible light
C-F-M	Compression forged method
CQ	Clinical Quality
BSS®	Balanced Salt Solution
SN	Serial Number

ABBREVIATION / SYMBOL	DESCRIPTION
	Unique Device Identification
GTIN	Global Trade Item Number
	Medical Device
	Catalogue number
	Batch code
	Use by: Year - Month
	Use by: Year - Month
	Manufacturer
	Country and date of manufacture
	Importer
	Distributor
	CE Mark with 4-digit number of the notified body
	Sterilised using irradiations
	Sterilised using steam
	Double sterile barrier system
	Single sterile barrier system
kGy	Kilogray
	Do not reuse
	Do not resterilize
	Consult instructions for use
	Consult electronic instructions for use
	Caution, consult accompanying documents
	Keep dry
	Keep away from sunlight
	Temperature limit (The applicable lower and upper temperature limits for storage are indicated on the device packaging. The applicable lower and upper temperature limits for transport are indicated on the outer packaging used for transportation.)

ABBREVIATION / SYMBOL	DESCRIPTION
	Do not use if packaging is damaged
	This side up
	Recyclable packaging
	Date of implantation
	Healthcare institution
	Patient information website
	Patient identification

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