PMMA Intraocular Lenses for Aphakia

INSTRUCTIONS FOR USE

Description

A PMMA intraocular lens (IOL) for Aphakia manufactured by OPHTEC is a high precision implant for replacing the human lens and is indicated for the surgical correction of the aphakic eye. The IOL consist of ultraviolet light absorbing polymethylmethacrylate (PMMA) with Tinuvin 326, a benzotriazole. Tinuvin 326 in an IOL exhibits effective UV filtration up to approximately 400 nm. Resolution of 20 D lens: > 70%. Light transmittance: 90%in the visible spectrum.

Indications

- Age-related cataract
- Traumatic cataract
- Congenital or juvenile cataract

Contraindications

- Only one eye with good vision
- Congenital bilateral cataract
- Recurrent or chronic iritis
- Rubella cataract
- Retina and optic nerve defects
- Corneal dystrophy (except in preparation for penetrating keratoplasty)
- Diabetic retinopathy
- Acute inflammation
- Severe iris atrophy
- Uncontrolled chronic glaucoma
- Vitreous loss or choroidal hemorrhage during surgery

Complications (lens related)

- Corneal endothelial dystrophy
- Hyphema
- Uncontrolled high intraocular pressure

Complications (not necessarily lens related)

- Intraocular infection
- Secondary cataract
- Corneal edema
- Corneal dystrophy
- Iris atrophy
- Pupillary block
- Pupillary membrane
- Iris prolaps
- Flat anterior chamber
- Retinal detachment
- Cystoid macular edema
- Iridocyclitis and vitritis
- Uveitis
- Glaucoma
- Hypopyon
- High ametropia and aniseikonia

Warnings

- Do not resterilize
- Do not reuse
- Do not use after the expiration date
- Do not use if the integrity of the product or the packaging has been compromised
- Do not store at temperatures > 40°C or 104°F
- Do not soak in fluids other than a sterile balanced salt solution

For single patient use only

Do not reuse, reprocess or resterilize, this may compromise the structural integrity of the device and/or create a risk of contamination. Malfunction of the device and cross-contamination may lead to injury, illness or death of the patient.

Precautions

- A high level of surgical skill is required for intraocular lens implantation
- For Artisan Aphakia: One or more iridotomies should be performed to reduce risk of pupillary block

Operating instructions

- Check the label for correct lens model, dioptric power and expiration date.
- Inspect the blister pack. Make sure it is not damaged.
- The IOL can become electrostatically charged and stick to the lid of the lens container. Tap slightly on the lid before opening the lens container.
- While keeping the container in a horizontal position, unscrew the cap and lift it.
- Grasp the lens gently with a forceps.
- Examine the lens carefully for damage or particulate matter under the microscope.
- Rinsing the IOL with sterile balanced salt solution will remove the electrostatic charge.

How supplied

The IOL is supplied sterile and dry in a lens container, which is sealed in a blister pack and placed in a box together with the identifying labels.

Returned lens policy

According to OPHTEC policy for the return/exchange/re-sterilization of intraocular lenses. Contact OPHTEC BV for more information.

Disclaimer of liability

OPHTEC BV shall not be liable for any injury or damage suffered by a patient as a result of:

- The surgical technique or implantation method used by the surgeon;
- Improper patient selection.

OPHTEC BV makes no expressed or implied warranties in connection with the resale or fitness for use of its intraocular lenses.

PRODUCT OF THE NETHERLANDS

Symbol	Explanation
	Batch Code
	Catalogue Number
	Sterilized by Ethylene Oxide
	Do Not Resterilize
	Do Not Reuse
	Use by (YYYY-MM: Year-Month)
	See Instructions for Use
	Serial Number
	Upper Limit of Temperature
	Do Not Use if Package is Damaged
	Manufacturer

 **OPHTEC BV**
Schweizerlaan 15
9728 NR Groningen
PO Box 398
9700 AJ Groningen
The Netherlands
Tel. +31 50 5251944
Fax +31 50 5254386
E-mail: info@ophtec.com

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