<i>CrownVisc® / Sodyum Hyaluronat</i>	TR
--	-----------

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User information:

CrownVisc 1.0% / 1.4% / 1.6% / 1.8% / 3.0% sodium hyaluronate. Viscoelastic solution in a single use glass syringe for intraocular use, and sterile cannula.

Composition Each ml contains:	CrownVisc 1.0%	CrownVisc 1.4%	CrownVisc 1.6%	CrownVisc 1.8%	CrownVisc 3.0%
Sodium hyaluronate	10,000 mg	14,000 mg	16,000 mg	18,000 mg	30,000 mg
Sodium chloride	8,500 mg	8,500 mg	8,500 mg	8,500 mg	8,500 mg
Disodium hydrogenphosphate	0,563 mg	0,563 mg	0,563 mg	0,563 mg	0,563 mg
Sodium dihydrogenphosphate	0,045 mg	0,045 mg	0,045 mg	0,045 mg	0,045 mg
Water for injection	q.s	q.s	q.s	q.s	q.s

Definition :

CrownVisc is a clear, colorless, odorless, sterile, isotonic, viscoelastic solution for intraocular use. CrownVisc is steam-sterilised under pressure in its sterile packaging. Cannula is sterilized using ethylene oxide. With its special viscoelastic properties, CrownVisc maintains a deep anterior chamber and protects the corneal endothelium throughout surgery. The viscoelastic properties of the preparation allow tissue structures to be moved without causing trauma, for example in cases of prolapse of the iris. Diagnostic and therapeutic procedures may be carried out in the eye without any impairment of vision. CrownVisc protects the intraocular tissue throughout surgery and prevents adhesion and synechia formation when performing ocular surgery. Once surgery is finished, CrownVisc must be removed completely by irrigation and/or aspiration. Any residue that has not been removed from the anterior chamber will be flushed out naturally via the trabecular system and Schlemm’s canal. If it is of not completely removed, intraocular pressure increases due to obstruction of the trabecular system.

Indications:

CrownVisc is indicated for use as an surgical aid in the anterior segment during cataract extraction, intraocular lens implantation and glaucoma surgery. CrownVisc maintains the anatomical space of the anterior segment with reduced trauma to the corneal endothelium during surgery. It is intended for intraocular use only.

Method of Administration and Posology:

CrownVisc is administered by a narrow-gauge blunt cannula. The amount required depends on the type of surgical procedure.

Contraindications:

CrownVisc must not be administered to patients who are known to be allergic to one of the constituents.

Interaction with other agents:

There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds such as benzalconium chloride solutions. CrownVisc should therefore not come into contact with operating instruments rinsed with these solutions or with Ophthalmic preparations containing ammonium compounds as a preservative. CrownVisc solution may appear cloudy or form precipitate when come into contact with quaternary ammonium compounds. When it is observed, the physician should remove the cloudy or precipitated material by irrigation and/or aspiration.

Warnings and Precautions:

This product should only be used by qualified people who are familiar with the use of such a product. CrownVisc is intended only for intra-ocular surgery. CrownVisc must not be resterilised. CrownVisc is intended for single use only. Do not administer after the expiry date. Opened and/ or damaged sterile packs must not be used. To prevent the formation of air bubbles, remove the cap from the syringe by twisting it off carefully. Care should be taken to prevent the introduction of trapped air bubbles into the anterior chamber.

Instructions for the correct removal of the cap:

Hold the Luer Lock Adapter as shown in ❶. Twist the cap carefully with the other hand in an anti-clockwise direction ❷. Then remove the cap as shown in ❸. Do not use a syringe if the cap is not on properly in its sterile packaging.

Instructions for the correct insertion of the cannula:

Hold the syringe as shown in ❹ Insert the enclosed 23-27 G cannula firmly as shown in ❺ (do not use any other cannula). Hold the cannula and lock it into position by twisting slightly in a clock-wise direction ❻. Hold CrownVisc during administration as shown in Figure 1. Avoid injecting an excessive volume of CrownVisc into the eye. There is a risk of the product drying out if it is left open for a long time during use. The product must be removed completely by irrigation and/or aspiration at the end of the operation. Otherwise, a mechanical blockage in the trabecular system could occasionally develop and cause a temporary postoperative rise in intraocular pressure. In patients diagnosed preoperatively with wide-angle glaucoma, pronounced myopia, diabetic retinopathy or uveitis, there is a high risk of an increase in intraocular pressure. If intraocular pressure rises to an excessive level, suitable treatment should be given to reduce the pressure. After use the cannula must be disposed of in a sharps container.

Storage:

Store CrownVisc at 2-25 °C (36-77 °F) in a cool, dry place, and protect from light, heat and frost.

Adverse Reactions: Unexpected adverse reactions and complications related to CrownVisc must be reported to ATAKAN DEDE - MİRAY MEDİKAL.

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