NEOCROWN® Cohesive / Sodium Hyaluronat TR

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Revizyon Numarası/Tarih: 04/23.01.2018

NEOCROWN® Cohesive / Sodium Hyaluronate EN

NEOCROWN Cohesive – Instructions for use

Product:

NEOCROWN, 1.4% sodium hyaluronate.

1,0 ml viscoelastic solution in a single-use glass syringe for intraocular use and one sterile cannula.

Composition:

1 ml solution contains 14,0 mg sodium hyaluronate, sodium chloride, disodium hydrogen phosphate, sodium dihydrogen phosphate and water for injection.

Description:

NEOCROWN is a clear, colourless and odourless, sterile, isotonic, viscoelastic solution for intraocular use. NEOCROWN is steam-sterilised under pressure in its sterile packaging. Cannula is sterilized using ethylene oxide. NEOCROWN contains sodium hyaluronate in a physiological saline solution with a pH of 6.0-8.0. The sodium hyaluronate contained in NEOCROWN is obtained from bacteria by fermentation, and is extremely well tolerated when used in the eyes. NEOCROWN is easy to instil. During injection, the viscosity decreases considerably in the cannula. After instillation, the original high viscosity is restored. With its special viscoelastic properties, NEOCROWN maintains a constant deep anterior chamber and protects the corneal endothelium throughout surgery. The viscoelastic properties of the preparation allow tissue structures to be moved without causing trauma, for example in cases of prolapse of the iris. Diagnostic and therapeutic procedures may be carried out in the eye without any impairment of vision. NEOCROWN protects the intraocular tissue throughout surgery and prevents adhesion and synechia formation when performing ocular surgery. At the end of the operation, NEOCROWN must be removed completely by irrigation and/or aspiration. Any residue that has not been removed from the anterior chamber will be flushed out naturally via the trabecular system and Schlemm's canal. If it is not completely removed, intraocular pressure increases due to obstruction of the trabecular system.

Indications:

NEOCROWN is indicated for use as an surgical aid in the anterior segment during cataract extraction, intraocular lens implantation and glaucoma surgery. NEOCROWN maintains the anatomical space of the anterior segment with reduced trauma to the corneal endothelium during surgery. It is intended for intraocular use only.

Administration and dosage:

NEOCROWN is administered through a narrow-gauge blunt cannula. The amount required depends on the type of surgical procedure.

Contraindications:

NEOCROWN must not be administered to patients known to be hypersensitive to any of the listed ingredients.

Interaction with other agents:

There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds such as benzalkonium chloride solutions. NEOCROWN should therefore not come into contact with surgical instruments rinsed with these solutions or with ophthalmic preparations containing ammonium compounds as a preservative. NEOCROWN solution may appear cloudy or form precipitate when come into contact with quaternary ammonium compounds. When it is observed, the physician should remove the cloudy or precipitated material by irrigation and/ or aspiration.

Warnings and Precautions:

This product may only be used by surgeons who have the requisite specialist knowledge and skill to handle it. NEOCROWN is intended only for intra-ocular surgery. NEOCROWN must not be re-sterilised. NEOCROWN is intended for single use only. Do not administer after the expiry date. Any sterile pack that has been opened and/or damaged must not be used. In order to prevent the formation of air bubbles, remove the cap from the syringe by twisting it off carefully. Care should be taken to prevent the introduction of trapped air bubbles into the anterior chamber.

Instructions for the correct removal of the cap:

Hold the Luer Lock Adapter as shown in **❶**. Twist the cap carefully with the other hand in an anti-clockwise direction **❷**. Then remove the cap as shown in **❸**. Do not use a syringe if the cap is open or is not on properly in its sterile packaging.

Instructions for the correct insertion of the cannula:

Hold the syringe as shown in **❶**. Insert the enclosed 27 G cannula firmly as shown in **❷** (do not use any other cannula!). Hold the cannula firmly and lock it into position by twisting slightly in a clock-wise direction **❸**. Hold NEOCROWN during administration as shown in Figure 1. Do not overfill the eye with NEOCROWN. There is a risk of the product drying out if it is left open during use over a long period of time. The product must be removed completely by irrigation and/or aspiration at the end of the operation. Otherwise, in some cases a mechanical blockage in the trabecular system could develop and cause a temporary post-operative rise in intraocular pressure. There is an elevated risk of an increase in intraocular pressure for patients who have been diagnosed prior to surgery with wide-angle glaucoma, severe myopia, diabetic retinopathy or uveitis. If intraocular pressure rises to an excessive level, appropriate treatment should be initiated to reduce it. After use, the cannula must be disposed of in a sharps container.

Storage:

Store NEOCROWN at 2 -8 °C (36°F -46 °F) in a cool, dry place, and protect from light, heat and frost.

Adverse Reactions: Unexpected adverse reactions and complications related to NEOCROWN must be reported to ATAKAN DEDE - MİRAY MEDİKAL.

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