Benutzerhandbuch / *User Manual*

MELAcontrol® Helix

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General guidelines

EN

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

| Symbols on the product | Description |
|---|---|
|  | Identifies a medical device |
|  | In affixing the CE mark, the manufacturer declares that this product fulfils the corresponding EU requirements. |

Intended use

The MELAcontrol Helix steam penetration test device is a hollow load process challenge device (PCD). The MELAcontrol Helix is designed for testing the air removal (steam penetration) capability of small steam sterilizers with type B cycles. This steam penetration test device may also be used as hollow devices in accordance with EN 285 for large steam sterilizers. The chemical indicator located in the capsule will show a defined colour change from yellow to dark blue/purple when exposed to a specific combination of time, steam and temperature. When all the air has been removed from the sterilization chamber, steam will penetrate through the steam penetration test device and the indicator will show a uniform colour change. If air or non-condensable gases are present the indicator will show distinct yellow markings. Each steam penetration test device is supplied with its own specific quantity of indicator strips.

PLEASE NOTE

MELAcontrol Helix supplements the physical measurements and is suitable for testing the following steam sterilizers:

- Steam sterilizers with type “B” cycles
- Steam sterilizers with type “S” cycles, of which the scope of supply covers the treatment of hollow bodies
- Large steam sterilizers in accordance with EN 285

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

Application

1. Unscrew the test body sealing cap.

2. Fold the indicator strip in the middle.

3. Insert the middle folded part of the indicator strip first into the sealing cap of the test body.

4. Guide the indicator strip into the sealing cap as far as possible.

5. Screw the sealing cap up to the stop in the test body.

CAUTION

Danger of burns from hot condensate! Position the test body for sterilization in a level position with the flat side facing downwards, otherwise hot condensate can run out of the test body after removal.

- Let the test body cool before removal.

PLEASE NOTE

The test body must be dry before every re-use; if possible blow it out with compressed air. Residue condensate in the test body has a negative effect on air removal.

- Insert an unused indicator strip in the plastic bracket of the sealing cap before every test.

Application for daily routine check

1. Place a new indicator strip in the test body.
2. Slide the test body onto a tray in the empty sterilization chamber.
3. Start the steam sterilizer in a *validated* sterilization program, according to the intended use.
4. Remove the test body from the steam sterilizer when the program run has ended.
5. Remove the indicator strip from the sealing cap. Evaluate the indicator strip.

Application as batch control

1. Place a new indicator strip in the test body.
2. Insert the test body in the sterilization chamber, together with the sterilization material.
3. Start the steam sterilizer in a *validated* sterilization program, according to the intended use.
4. Remove the test body from the steam sterilizer when the program run has ended.
5. Remove the indicator strip from the sealing cap. Evaluate the indicator strip.

PLEASE NOTE

Depending on the manufacturing processes, the colour change can be different (dark blue or light violet). The decisive factor for the evaluation is the consistent and complete colour change of the indicator strip whereas a slight gradient to the middle is still acceptable.

Evaluation

| | | |
|--|--|---|
| | completely changed indicator strip | Sterilization successful sufficient air removal and steam penetration |
| | changed indicator strip with slight gradient to the middle | Sterilization successful sufficient air removal and steam penetration |
| | incompletely changed indicator strip | Sterilization not successful insufficient air removal and steam penetration |

Documentation

The evaluation of the results must be performed by a sufficiently trained person. The colour change is caused by an irreversible chemical reaction. If stored correctly, the indicator strips can be archived over years. EN ISO 15882 only requires a description of the results. This can be performed in written or digital form. Comply with the valid national regulations.

Storage

The indicators must be stored in a cool, dark and dry place before and after use. Do not store the indicator strips together with chemical substances and gases (e.g. ethylene oxide [EO]), hydrogen peroxide [H₂O₂] or instruments sterilized with H₂O₂. Incorrect storage can cause undesirable colour changes. Further information regarding storage is provided on the packaging. A secure determination of failure in the process can only be guaranteed up to the specified expiring date.

Important information regarding routine operation

After 250 indicator have been used the product must be replaced.

Please comply with the current recommendations issued by the Robert Koch Institute (RKI) and the information contained in DIN 58946-7. Further information regarding routine operation is to be found in the user manual of your steam sterilizer.