<i>CrownGel® / Hidroksipropil Metilselüloz</i>	TR
---	-----------

Ürün:

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Product:

CrownGel 2.0% - 2.4% - Instructions For Use

CrownGel, 2.0 ml viscoelastic solution in a single use glass syringe for intraocular use, and sterile cannula.

Composition 1mL solution contains	CrownGel 2.0%	CrownGel 2.4%
Hydroxypropylmethylcellulose	20.00 mg	24.00 mg
Sodium Chloride	6.40 mg	6.40 mg
Potassium Chloride	0.75 mg	0.75 mg
Calcium Chloride 2H ₂ O	0.48 mg	0.48 mg
Magnesium Chloride 6H ₂ O	0.30 mg	0.30 mg
Sodium Acetate 3H ₂ O	3.90 mg	3.90 mg
Sodium Citrate 2H ₂ O	1.70 mg	1.70 mg
Other Components	Sufficient quantity of HCl or NaOH, Water for injection	

Description:

CrownGel is an isotonic, clear, viscoelastic and sterile solution for intraocular application. The osmolality of CrownGel is 270 to 400 mOsm/kg; the pH is 6.8 to 7.6. The CrownGel has an average molecular weight of 480000 g/mol. CrownGel is steam-sterilised under pressure in its sterile packaging. Cannula is sterilized using ethylene oxide. As an auxiliary device during surgery of the anterior segment of the eye, CrownGel maintains the depth of the anterior chamber and protects the surrounding intraocular tissue. CrownGel contains 2.0% or 2.4% hydroxypropyl methylcellulose and supports the intraocular manipulations due to its viscoelastic and moistening properties. CrownGel does not cause any inflammatory or immunogenic reactions. CrownGel is free of preservatives. Due to its solubility in water, CrownGel can be removed at the end of surgery simply by irrigation and aspiration.

Indications:

CrownGel serves as a volume substitute for aqueous humor during intraocular operations, such as extraction and insertion of an intraocular lens.

CrownGel maintains the depth of the anterior chamber during surgery and reduces the danger of traumatization of the corneal endothelium, the iris and ciliary body by direct contact with surgical instruments.

Method of Administration and Posology:

CrownGel is carefully injected into the anterior chamber using the supplied sterile single-use 23G cannula. It is recommended for intraocular lens implantation to cover the implant and instruments with CrownGel immediately before surgery. This measure is to additionally protect the endothelium and the surrounding tissue. CrownGel volume to be injected depends on the type of surgery. As a substitute for viscoelastic lost by flow or irrigation, CrownGel can be injected several times. At the end of surgey, CrownGel must be removed completely using a suitable irrigation/aspiration device.

Contraindications:

CrownGel should not be used in cases of a hypersensitivity against hydroxypropyl methylcellulose.

Warning and Precautions:

CrownGel is designed for use by ophthalmologists only. All normal precautions for intraocular surgery must be observed. CrownGel and the supplied sterile single-use cannula are designed for single intraocular application only. Do not use CrownGel after the expiration date. Any sterile pack that has been opened and/or damaged must not be used. In order to prevent the formation of air bubbles, remove the cap from the syringe by twisting it off carefully. Prior to removing the cannula from the eye, do not exert any more pressure into the plunger. In this way aspiration of an air bubble into the cannula can be avoided. Avoid injecting an excessive volume of CrownGel into the eye. There is a risk of the product drying out if it is left open for a long time during use. The product must be removed completely by irrigation and/or aspiration at the end of the operation. Otherwise, a mechanical blockage in the trabecular system could occasionally develop and cause a temporary post-operative rise in intraocular pressure. In patients diagnosed preoperatively with wide-angle glaucoma, severe myopia, diabetic retinopathy or uveitis, there is a high risk of an increase in intraocular pressure. If intraocular pressure rises to an excessive level, suitable treatment should be given to reduce the pressure. After use, dispose of the product with medical waste in accordance with local waste management regulations.

Instructions for the correct removal of the cap:

Hold the Luer Lock Adapter as shown in ❶. Twist the cap carefully with the other hand in an anti-clockwise direction ❷. Then remove the cap as shown in ❸. Do not use the syringe if the cap is not on properly in its sterile packaging.

Instructions for the correct insertion of the cannula:

Hold the syringe as shown in ❹. Insert the enclosed 23 G cannula firmly as shown in ❺ (do not use any other cannula). Hold the cannula and lock it into position by twisting slightly in a clock-wise direction ❻. Hold CrownGel during administration as shown in Figure 1. After use the cannula must be disposed of in a sharps container.

Storage:

Store CrownGel at 2 -25 °C (36°F -77 °F) in a cool, dry place, and protect from light, heat and frost.

Adverse Reactions:

Unexpected side effects and complications related to product should be reported to the competent authority of the member state and to ATAKAN DEDE-MİRAY MEDİKAL.

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