[REF] 70002 - 70009

EN	ENGLISH
----	---------

PVA Eye Spears

DESCRIPTION

S_{IDA}-SPEARS PVA Eye Spear is a sterile, single use device consisted of a spear and a handle. The spear is made of PVA (polyvinyl alcohol) material which can absorb aqueous fluids several times of its mass instantaneously upon contact with the fluids. The handle is made of polypropylene plastic material.

INDICATIONS –INTENDED USE

The PVA Eye Spear is indicated for use in ophthalmic or other microsurgical procedures for absorbing or mopping fluids from the operative field.

INTENDED PATIENT POPULATION

Male and female older than 7 years.

INTENDED USER

Healthcare professionals (mainly eye surgeons) on the field of ophthalmology.

CONTRAINDICATIONS

- Absolute Contraindication: Patient allergy to PVA (polyvinyl alcohol) materials.
- Relative Contraindication: Pregnant women or patients who have had any adverse reactions by using PVA related materials.

SIDE EFFECTS

None.

PRECAUTIONS

- Sterility is guaranteed unless unit package is opened or damaged.
- S_{IDA}-SPEARS should not be moistened with an alcoholic solution, as then the structure of the material will change.

USE PROCEDURE

The product is packed inside an aluminum foiled pouch which in turn is sealed in a pouch or plastic bag (depending on the model). The pouch consists of a Tyvek (synthetic paper) top layer and a transparent plastic bottom layer.

- To open the outer packaging, gently pull open to separate the two layers. In case of bag, hold it upright, open it by hands or carefully by using a pair of scissors.
- Empty the foiled pouch onto a sterile field.
Caution: Maintain sterile technique when contacting the inner foiled pouches.
- While maintaining sterile technique, tear open the foiled pouch.
- Remove spears and place them in a sterile field.
Caution: In order to avoid inadvertent contact of the spear with aqueous fluids, keep the spears in a dry place within the sterile field.
- By holding the plastic handle, gently touch the tip of the spear at the site where excess fluid is to be removed.
- Hold still for approximately 3 seconds to allow for complete swollen of the PVA material.
Caution: Do not squeeze or press the PVA eye spear against the surface of the surgical sites.
- Dispose of the used spear into a biological waste container.
- When necessary, use additional spears to absorb fluid until all excess fluid is satisfactorily removed.

WARNINGS

- For single use only. Do not reuse, as it may result in compromised device performance and risk of cross contamination.
- Do not attempt to use the spear for cleaning of the surgical sites or removing debris, such as vitreous material that might be generated during the surgery.
- The Eye Spear should not be used as an applicator for applying medication to the surgical sites. Also, it shall not be used as a hemostasis tool for compressing the spear against the bleeding site.

PRESERVATION & HANDLING

- The product should be preserved in a dry, well-ventilated place and in sealed conditions.
- Do not use if package is opened or damaged. The identification tag with the Lot no. should not be separated while in preservation.
- Storage conditions: 0°C↔ 37°C.
- Dispose of used sponges according to local medical regulations.

REPORTING

Adverse events and/or complications that may reasonably be regarded as product-related should be reported to the manufacturer.

DE	DEUTSCH
----	---------

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