
RELEASED TO
PRODUCTION

Qube®
Qube® pro
PentaSys
PentaSys2

Bi-Blade®

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CE 0197

DE

Directions for Use Sterile Bi-Blade® Vitreous Cutter

The Bi-Blade® vitreous cutter is used to remove vitreous and intraocular tissue. The Bi-Blade® vitreous cutter is labeled for use with the RUCK ophthalmic surgical systems Qube®, Qube® pro, PentaSys and PentaSys2.

STERILE R

⚠️ WARNUNG

Sterility not guaranteed if package has been opened or damaged.
This product is for single patient use only.
Do not re-sterilize.
Risk of cross-contamination if product is re-used.

⚠️ CAUTION

Federal (USA) law restricts this device to use or sale by or on order of a physician.

⚠️ CAUTION

344 KPa (50 psi) Maximum Pressure at Tubing End

⚠️ WARNUNG

Test Vitreous Cutter before use.
If bubbles are seen at cutting port, do not use cutter.

⚠️ CAUTION

Do not operate the cutter dry since the cutting action is dependent on a liquid environment.

⚠️ CAUTION

Failure to properly attach the air supply tubing to the appropriate pressure source will affect the operation of the vitrectomy cutter.

⚠️ CAUTION

Aspiration flow may be greater than expected. Start with a reduced vacuum level.

INSTALLATION

- Tubing set is color coded for easy identification; GREEN striped line for cutting actuation pressure and CLEAR tubing for aspiration.
- Handle using appropriate aseptic technique.
- Attach GREEN striped tubing to the vitrectomy actuation connector on system console.
- Attach CLEAR tubing to the aspiration system of vitrectomy system console.

TESTING AND PRIMING

Note: Refer to user manual of system console for additional instruction on priming and testing the cutter.

- Place needle of cutter in a container of sterile solution and actuate cutting. Look for bubbles or other leaks. If bubbles are seen at cutting port, do not use cutter.
- Under the microscope, test for appropriate cutting action of the vitrectomy cutter. Do not run the cutter outside sterile solution for more than a few moments. If the vitrectomy cutter is not operating properly, do not use cutter.
- Place needle of cutter in a container of sterile solution and prime tubing set. To prime, actuate aspiration until all bubbles are removed from tubing set.

Glossary of Symbols

Symbol	Title/Description	Applicable Standard	Reference Number
	Manufacturer	ISO 15223-1 BS EN 980	5.1.1 5.12
	Authorized Representative in the European Community	ISO 15223-1 BS EN 980	5.1.2 5.13
	Use-by date	ISO 15223-1 BS EN 980	5.1.4 5.3
	Batch code	ISO 15223-1 BS EN 980	5.1.5 5.4
	Catalogue number	ISO 15223-1 BS EN 980	5.1.6 5.10
	Sterilized using irradiation	ISO 15223-1 BS EN 980	5.2.4 5.8.3
	Do not re-sterilize	ISO 15223-1 BS EN 980	5.2.6 5.22
	Do not re-use	ISO 15223-1 BS EN 980	5.4.2 5.2
	Do not use if package is damaged	ISO 15223-1 BS EN 980	5.2.8 6.3
	Consult instructions for use	ISO 15223-1 BS EN 980	5.4.3 5.18
	Caution	ISO 15223-1 BS EN 980	5.4.4 5.11
	Caution: Federal (U.S.) Law restricts this device to sale by or on the order of a physician.	21 CFR 801	§ 801.109(b)(1)
	European Conformity – CE Mark	Council Directive 93/42/EEC	Article 17
	Contains or presence of phthalates (DEHP)	EN 15986	4.2

If any MID Labs products are received in a defective condition, your distributor or MID Labs is to be notified immediately.

Medical Instrument Development Laboratories, Inc.
557 McCormick Street
San Leandro, CA 94577 U.S.A.
Phone: 1-800-929-5227 (Inside U.S. only)
Phone: 510-357-3952
Fax: 510-357-1582
Email: info@midlabs.com

Each MID Labs product is identified by a lot number. Refer to this number when discussing the product with your MID Labs representative.

Use of accessories other than those under the branding of Qline® may affect system performance and create potential hazards.

Manufacturer assumes no responsibility for complications that may arise from use of accessories other than those under Qline® branding or improper cleaning, sterilization, or use, of Qline® products.

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