<i>Ocudeka Decaline</i>	TR
-------------------------	----

Açıklama

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Description

Ocudeka Decaline is a high density bicyclical fluorocarbon compoud that is composed only of C-C and C-F bonds.

Pyhsical Prosperities

Ocudeka Decaline	
Formula	: C10F18
Molecular weight	: 462,08
Density	: 1,908 g / cm³ at 20 °C
Purity	: 100 % fluorinated perfluorocarbon
Boiling point	: 142 °C

Instructions For Use

Ocudeka Decaline is per-operative devices for vitreoretinal surgey. Following a partial or complete vitrectomy, Ocudeka Decaline liquid is injected exactly opposite the optical disc. Due to the lower viscosity standard cannula from 20 to 23 G can be used. In the case of large tears penetration of the compound beneath the retinal tissue must be avoided. At the end of operation Ocudeka Decaline must be completely removed from the eye.

Indications

As a agent for:

- Retinal detachment
- Proliferative diabetic retinophy
- Macular holes with peripheral RD
- Retrieval of foreign bodies from vitreous
- Giant tears
- Trauma
- Laser coagulation and cryotherapy
- Refloating luxated lenses
- Short-term stabilising of reattached retinas

Contraindications

Due to the high density it cannot be excluded that Ocudeka Decaline may induce changes in the retina or choroid with long-term use Ocudeka Decaline should therefore only be used intraoperatively and for a short time.

Side Effects

Side effects are not know.

Presentation And Sterilisation

Ocudeka Decaline can be supplied in a sterile and pyrogen-free vial of 5ml, 7ml and 10ml. Vials are steam sterilized. Supplied with one 5 ml sterile disposable syringe. Disposable syringe is sterilized using ethylene oxide.

Storage Conditions

Ocudeka Decaline must be stored between 2°C - 30°C (36-86°F).

Warnings

Ocudeka Decaline is designed for use by ophthalmologists only. Do not use this product after expiry date. Only use undamaged packaging. Vials are only intended for single use on one patient.

Adverse Reactions

Unexpected adverse reactions and complications related to Ocudeka Decaline must be reported to ATAKAN DEDE – MİRAY MEDİKAL.

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