GOT Filtration and injection

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GOT Filtration and injection

WARNING: This instruction sheet contains important information and warnings on the correct use of the medical device. Carefully read it before starting any operation and closely follow the directions for use.

COMPOSITION AND PACKAGING

GOT Filtration and injection is composed of a PTFE 0.22 µm filter with a PVC connector. Each box includes 30G needles and patient wristbands.

INDICATIONS FOR USE

GOT Filtration and injection, accessory to filter and inject ophthalmic gas, is intended for the exclusive use together with the ophthalmic gases SF₆, C₂F₆, C₃F₈, contained in systems provided with 4 mm valve spindle and/or male luer lock connectors that operates at pressure up to 1 bar, and sterile syringes with luer connector.

GOT Filtration and injection is intended for use at centres where vitreoretinal surgeries are performed. Consequently, its use is restricted to qualified medical personnel only. Any other use must be considered improper.

DIRECTIONS FOR USE

- Make sure that the package is intact.
- Remove the **GOT Filtration and injection** from the pouch in the sterile field.
- Connect the filter (2) to a sterile syringe with a luer connector under aseptic conditions.
- Attach the connector (1) to the cylinder containing the ophthalmic gas, according to instructions supplied by the gas manufacturer.
- Take the gas according to the instructions supplied by the gas manufacturer.
- Disconnect the filter (2) and the connector (1) from the syringe.
- Immediately mount the needle to the syringe.

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- The product is ready for the intraocular injection.

WARNINGS AND PRECAUTIONS

- Before use, store as indicated on label.
- Integrity is warranted only if package is intact – Do not use if package is open or damaged or any component is received in a defective condition.
- The expiry date is shown on product label – Do not use after the expiry date.
- Single use only – Do not re-sterilize. Use or sterilization of single use products that have been already used may cause cross-contamination.
- Product shall be used by qualified medical personnel only, in surgery room or in any case in the sterile field.
- The device is intended for the exclusive use together with the ophthalmic gases SF₆, C₂F₆, C₃F₈, contained in systems provided with 4 mm valve and/or male luer lock that operates at pressure up to 1 bar, and sterile syringes with luer connector.
- Absolutely do not remove the filter from the syringe during gas withdrawal and possible mixture with air.

- Patient implanted with intraocular gas must be invited to wear for an appropriate period of time the enclosed wristband containing important restrictions, including limitations on the use of nitrous oxide (N20) in subsequent surgical or dental procedures and exposition to pressure variations (elevation, travel by plane, diving with or without cylinders, hyperbaric chamber treatment).**

- After use, components in direct contact with patient body fluids are potentially bio-hazardous – Take the necessary steps for self-protection and dispose of in compliance with national regulations on workplace safety and quality assurance principles.
- Dispose of other components and empty packages in compliance with regulations in force – Avoid PVC dispersion into the environment.

STERILIZATION

GOT Filtration and injection and needle are sterile by ethylene oxide.

Single use only – Do not resterilize.

STORAGE

Store at room temperature.

For any report on this product, quote the lot number on package. For any comments or suggestions, go to www.alchimiasrl.com.

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