<i>Sodyum Hyaluronat OVD</i>	TR
-------------------------------------	-----------

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Usage information:

CrownVisc/NanoVisc 1.0% / 1.4% / 1.6% / 1.8% / 3.0% sodium hyaluronate. Viscoelastic solution in a single use glass syringe for intraocular use, and a sterile cannula.

Composition Each ml contains:	1.0%	1.4%	1.6%	1.8%	3.0%
Sodium hyaluronate	10 mg	14 mg	16 mg	18 mg	30 mg
Sodium chloride	8,5 mg	8,5 mg	8,5 mg	8,5 mg	8,5 mg
Disodium hydrogenphosphate	0,563 mg	0,563 mg	0,563 mg	0,563 mg	0,563 mg
Sodium dihydrogenphosphate	0,045 mg	0,045 mg	0,045 mg	0,045 mg	0,045 mg
Water for injection	q.s	q.s	q.s	q.s	q.s
Viscosity Ranges	200000 -400000 mPa.s	300000 - 470000 mPa.s	350000 - 550000 mPa.s	450000 - 750000 mPa.s	350000 - 650000 mPa.s
Molecular Weights	3.8 - 4.5 MDa	2.9 - 3.8 MDa	2.9 - 3.8 MDa	2.9 - 3.8 MDa	1.8 - 2.5 MDa
pH	6.8 – 7.6	6.8 – 7.6	6.8 – 7.6	6.8 – 7.6	6.8 – 7.6
Osmolality	270-400 mOsm/kg	270-400 mOsm/kg	270-400 mOsm/kg	270-400 mOsm/kg	270-400 mOsm/kg

Definition :

CrownVisc/NanoVisc is a clear, colorless, odorless, sterile and isotonic viscoelastic solution for intraocular use. CrownVisc/NanoVisc is steam-sterilised under pressure in its sterile packaging. Thanks to its special viscoelastic properties, CrownVisc/NanoVisc maintains a deep anterior chamber and protects the corneal endothelium throughout the surgical operation. The viscoelastic properties of the preparation allow tissue structures to be moved without causing trauma, for example in cases of prolapse of the iris. Diagnostic and therapeutic procedures may be carried out in the eye without any impairment of vision. CrownVisc/NanoVisc protects the intraocular tissue throughout surgery and prevents adhesion and synechia formation when performing ocular surgery. Once surgery is finished, CrownVisc/NanoVisc must be removed completely by irrigation and/or aspiration. Any residue that has not been removed from the anterior chamber will be flushed out naturally via the trabecular system and Schlemm’s canal. If it is not completely removed, intraocular pressure may increase due to obstruction of the trabecular system. CrownVisc/NanoVisc 1.0%, 1.4%, 1.6%, 1.8% shows cohesive properties, CrownVisc/NanoVisc 3.0% shows dispersive properties.

Cannula: Sterimedix 25g x 7/8" (0.5 x 22mm) / 27g x 7/8" (0.4 x 22mm) The visco-elastic cannula is produced for single use for ophthalmic irrigation/aspiration procedures. The cannula is sterilized with ethylene oxide.

Indications:

CrownVisc/NanoVisc is indicated for use as a surgical aid in the anterior segment during various intraocular surgical procedures including cataract extraction, intraocular lens implantation and glaucoma surgery. CrownVisc/NanoVisc maintains the anatomical space of the anterior chamber with reduced trauma to the corneal endothelium during surgery. It is intended for intraocular use only.

Method of Administration and Posology:

CrownVisc/NanoVisc is administered by a narrow-gauge blunt cannula. The amount required depends on the type of surgical procedure.

Contraindications:

Care should be taken in patients with hypersensitivity to any components of CrownVisc/NanoVisc.

Interaction with other agents:

There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds such as benzalkonium chloride solutions. CrownVisc/NanoVisc should therefore not come into contact with surgical instruments rinsed with these solutions or with ophthalmic preparations containing quaternary ammonium compounds as a preservative. CrownVisc/NanoVisc solution may appear cloudy or form precipitate when come into contact with quaternary ammonium compounds. When it is observed, the physician should remove the cloudy or precipitated material by irrigation and/or aspiration.

Warnings and Precautions:

CrownVisc/NanoVisc is designed for use by ophthalmologists only. CrownVisc/NanoVisc is intended only for intraocular surgery. CrownVisc/NanoVisc must not be resterilised. CrownVisc/NanoVisc is intended for single use only. Do not use after the expiry date. Opened and/or damaged sterile packs must not be used. To prevent the formation of air bubbles, remove the cap from the syringe by twisting it off carefully. Care should be taken to prevent the introduction of trapped air bubbles into the anterior chamber. Use only if the solution is clear. After use, dispose of the product with medical waste in accordance with local waste management regulations.

Instructions for the correct removal of the cap:

Hold the Luer Lock Adapter as shown in **❶**. Twist the cap carefully with the other hand in an anti-clockwise direction **❷**. Then remove the cap as shown in **❸**. Do not use the syringe if the cap is not on properly in its sterile packaging.

Instructions for the correct insertion of the cannula:

Hold the syringe as shown in **❹** Insert the supplied 25-27G cannula firmly as shown in **❺**. (do not use any other cannula). Hold the cannula and lock it into position by twisting slightly in a clock-wise direction **❻**. Hold CrownVisc/NanoVisc during administration as shown in Figure 1. Avoid injecting an excessive volume of CrownVisc/NanoVisc into the eye. There is a risk of the product drying out if it is left open for a long time during use. The product must be removed completely by irrigation and/or aspiration at the end of the operation. Otherwise, a mechanical blockage in the trabecular system could occasionally develop and cause a temporary postoperative rise in intraocular pressure. In patients diagnosed preoperatively with wide-angle glaucoma, pronounced myopia, diabetic retinopathy or uveitis, there is a high risk of an increase in intraocular pressure. If intraocular pressure rises to an excessive level, suitable treatment should be given to reduce the pressure. After use the cannula must be disposed of in a sharps container.

Storage:

Store CrownVisc/NanoVisc at 2-25 °C (36-77 °F) in a cool, dry place, and protect from light, heat and frost.

Adverse Reactions: Unexpected side effects and complications related to product should be reported to the competent authority of the member state and to ATAKAN DEDE-MİRAY MEDİKAL.

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