ES

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Synergy Plus 1.8%

Hialuronato de Sodio

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EN Instructions For Use

Synergy Plus 1.8%
Sodium hyaluronate

Produit:
Synergy Plus 1.8% is a viscoelastic solution containing sodium hyaluronate for injection during eye surgery.

Contents:
Each pack contains 1.1 ml of viscoelastic solution supplied in a single-use sterile glass syringe, for intraocular use, sealed in a blister pack.

Composition:
1.1 ml of Synergy Plus 1.8% contains sodium chloride, disodium hydrogen phosphate dodecahydrate, sodium dihydrogen phosphate dihydrate, Taurine, Sodium Hyaluronate (18mg/ml) and water for injection.

Description:
The Synergy Plus 1.8% viscosurgical device is a clear, sterile, non-pyrogenic buffered saline solution containing highly purified, and high-molecular-weight sodium hyaluronate obtained by bacterial biofermentation.
Synergy Plus 1.8% has a pH of 6.8 to 7.6 and osmolality of 300 to 350 mOsm/kg, similar to the aqueous humour.

Characteristics:
Synergy Plus 1.8% is a biocompatible solution with viscoelastic properties that lubricate, support and protect the tissues of the eye during ophthalmic surgery by forming a protective layer. Synergy Plus 1.8% maintains anterior chamber depth during surgery and facilitates procedures, with less trauma to the corneal endothelium and other surrounding tissues. Synergy Plus 1.8% helps create a clear visual field for inspection of the anterior segment and retina using contact lenses for examination. Synergy Plus 1.8% does not interfere with re-epithelialisation or normal wound healing.
Synergy Plus 1.8% is easy to inject. In fact, even though sodium hyaluronate is viscous at rest, its pseudoplastic properties ensure good flow when pressure is applied to inject the solution through the cannula. After the injection, the viscoelastic properties are immediately restored.
Synergy Plus 1.8% is free from preservatives and is easy to remove by adequate irrigation and/or aspiration on completion of all surgical procedures.

Indications:
Synergy Plus 1.8% is indicated in eye surgery of the anterior segment, including cataract removal, insertion of intraocular lenses, corneal transplant, glaucoma and retinal detachment. Synergy Plus 1.8% can be used to protect and lubricate intraocular lenses and injection systems before implantation.

Dosage and administration:
Cataract surgery – Intraocular lens insertion
Before extracting the lens, Synergy Plus 1.8% is introduced into the anterior chamber to protect the corneal endothelium and maintain a deep anterior chamber.

Corneal transplant surgery
On removal of the recipient's cornea, the anterior chamber is filled with Synergy Plus 1.8%. In transplantation, the donated corneal graft can be placed on the surface of the solution and sutured into place. Additional Synergy Plus 1.8% may be injected to maintain a deep anterior chamber. Synergy Plus 1.8% can also be used to protect the endothelial layer that is exposed while the graft is being prepared for the recipient.

Glaucoma surgery
Synergy Plus 1.8% is injected before a trabeculectomy to reconstruct the anterior chamber. On completion of the procedure, further injection of Synergy Plus 1.8% may be needed to achieve good subconjunctival filtration and prevent tissue adhesion.

Retinal detachment
Synergy Plus 1.8% should be slowly introduced into the vitreous cavity. With caution and by modulating the injection, Synergy Plus 1.8% can be used to separate the membranes from the retina for safe excision and release of traction. Synergy Plus 1.8% also helps move tissues into

the desired positions. For example, it is useful for gently pushing back a detached retina, unrolling a retinal flap or holding the retina against the sclera during repositioning.

Contraindications:
At present, there are no known contraindications associated with the use of Synergy Plus 1.8% when used as recommended. Care must be taken when using the product in patients who are sensitive to any of the components in this solution.

Warnings and precautions:
Any viscoelastic residue not removed at the end of the surgical procedure will be eliminated naturally via the trabecular system and Schlemm's canal. However, there is a risk of blockage of the drainage channels, which could lead to increased intraocular pressure. To avoid the risk of increased intraocular pressure, irrigation/aspiration is used to eliminate product residue on completion of surgery. In the event of increased intraocular pressure, treat with the most appropriate therapy.
Patients who have been diagnosed with open-angle glaucoma, severe myopia, diabetic retinopathy or uveitis are more likely to receive any treatment.

Use of a 27G sterile cannula is recommended.
Do not reuse the syringe or cannula. Reuse of these carries the risk of contamination and infection for the patient.
Do not resterilize the pre-filled syringe.
Do not use the same package for more than one patient.
Do not use if the package is damaged or open.
Do not use after the expiry date printed on the packaging.
Dispose of the syringe in accordance with established medical practice and applicable local requirements.
There is no evidence on the safety of Synergy Plus 1.8% during pregnancy and breastfeeding. Use during pregnancy and breastfeeding must be at the discretion of the surgeon.


Interactions:
Sodium hyaluronate is incompatible with quaternary ammonium compounds such as benzalkonium chloride solutions. Therefore, syringes pre-filled with Synergy Plus 1.8% should never come into contact with surgical instruments rinsed using these solutions or with ophthalmic products containing quaternary ammonium compounds as preservatives.

Instructions for use:
FOR INTRAOCULAR USE. FOR SINGLE USE ONLY.
Synergy Plus 1.8% should only be used by specialists trained in intraocular surgery. The amount of product to be used should be adapted to the type of surgical procedure.
It is recommended that a small amount of Synergy Plus 1.8% be drawn from the syringe and carefully inspected before using the product.
It is recommended that the syringe be carefully inspected before using the product.
It is recommended that excessive pressure on the plunger be avoided.


Sterile opening technique
Peel back the Tyevek covering of the blister.
Bend the blister towards the middle so that the plastic plunger of the syringe comes fully out of its housing.
Remove the syringe from the blister and place it in the sterile field.


Assembly
Remove the plastic cap from the syringe tip by unscrewing it anti-clockwise.

Attach the cannula to the syringe and ensure that it is securely inserted, screwing it on clockwise.

Verify the proper functioning of the device by gently depressing the plunger until Synergy Plus 1.8% comes out of the tip of the cannula.

Storage:
Store Synergy Plus 1.8% between 2°C and 25°C. Do not freeze. If stored in a refrigerator at low temperature, Synergy Plus 1.8% must be left to stand at room temperature for at least 20 minutes before use.
Protect from light and impact.

Symbols:	
	Manufacturer
	Do not reuse
	Do not resterilize
	Sterilized using steam
	Consult instructions for use
	Caution
	Keep away from sunlight
	Keep dry
	Do not use if package is damaged

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