LABTICIAN RETINAL IMPLANTS

Description

Labtician silicone retinal implants are manufactured using 100% implant grade silicone elastomers. All Labtician retinal implants are MRI compatible and are latex free. Labtician Retinal Implants (Exopplants) are designed to be used alone or in combinations which offer the best possible shapes and sizes to effectively complete the scleral buckling procedure.

How Supplied

Each Labtician Retinal Implant is supplied sterile, in a sealed, double pouch format, five units per carton. Labtician Retinal Implants are ethylene oxide sterilized. The external surfaces of the outer pouch are not sterile.

Indications

Retinal Implants are used in the surgical treatment of retinal detachments. The implant creates an indentation in the sclera and choroid which approximates the retina to the pigment epithelium. This procedure is accomplished by local indentation or through the use of an implant that encircles the globe. Circling Bands and Silicone Strips encircle the eye and are used alone or in combination with other Implants. Grooved Strips are designed for narrow or wide beds or for high scleral buckles. Tire shaped Implants and Wedges are used under encircling elements for wide scleral beds; for breaks near the ora serrata; for multiple breaks; and for high scleral buckles. Sleeves are used to secure the encircling elements. Buttons are indicated for use under grooved strips to produce higher buckles. Meridional Implants give additional buckling in the meridional direction, and are used under curved or grooved implants. The Pad is used to close scleral ruptures. Silicone Sponge Implants are used alone or in combination and are often indicated for segmental buckling procedures.

Contraindications

Labtician retinal implants should not be used for patients who have a known sensitivity to silicone polymers.

Adverse Reactions

No surgical implant material is completely free of adverse reactions in the human body. However, long term experience with Labtician Retinal Implants has shown that the risk of a reaction is very low. As with all retinal implants, the possibility of extrusion, erosion or inadvertent infection exists. Product to be used only by qualified surgeons.

Warnings

IGNORING THESE WARNINGS OR REUSING PRODUCT MAY RESULT IN NON-STERILE OR DAMAGED PRODUCT.

- Do not autoclave or attempt to re-sterilize by any method.
- Do not use if sterile package integrity has been compromised

Re-sterilization

If an implant is prepared for implantation but is not subsequently used, the implant should not be re-sterilized or otherwise used. It should be discarded.

NEVER ATTEMPT TO RESTERILIZE LABTICIAN RETINAL IMPLANTS.

Additional Risks

Although rare, risks that may be associated with implantation of scleral buckles include: Infection, Visual Distortion (shortened axial length, myopia), Strabismus (ocular motility or torsion), Retinal Ischemia, Vitreoretinal Incarceration, Choroidal Ischemia, Extrusion (external erosion), Transscleral (internal) Erosion, Hypotony (secondary to internal erosion), limitation of eye growth/development of amblyopia (in children), and procedure failure. The benefits of the device (reattachment of the retina, prevention of blindness) outweigh the residual risks.

Disposal

Disposal of damage or used device should be done as per user's country disposal regulations and hospital/clinic protocol.

Disclaimer of Liability

Labtician Ophthalmics, Inc. shall not be liable for any injury or damage suffered by a patient as a result of:

- Any implantation method or technique used by a physician to implant this retinal implant.
- Any selection and use of the retinal implant for any particular patient or patient's condition. Labtician Ophthalmics, Inc. makes no expressed or implied warranties in connection with the sale of its retinal implants and specifically disclaims warranty liability upon merchantability or fitness for use.

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© 2009
Labtician Ophthalmics, Inc. não será responsável por qualquer lesão ou dano sofrido por um paciente como resultado de: