IXIUM TWIN

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FRANÇAIS

DESCRIPTION :

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Notice IXIUM TWIN - M0033 - 2016-07.

PRODUCT DESCRIPTION :

IXIUM TWIN is a viscoelastic solution of high molecular weight sodium hyaluronate, highly purified, clear, isotonic, sterile and pyrogen-free, for intraocular injection in surgery of the anterior segment of the eye. **IXIUM TWIN** is supplied in Luer Lock 2.0-mL prefilled glass syringes with a single-use disposable injection cannula.

CHARACTERISTICS :

The sodium hyaluronate used in the manufacture of **IXIUM TWIN** is a pharmaceutical-grade high molecular weight polysaccharide consisting of sodium glucuronate and N-acetylglucosamine obtained by genetic engineering (bacterial fermentation). Sodium hyaluronate is a physiological substance present in large quantities in numerous connective tissues in man and animals, and, particularly, in the vitreous humor, synovial fluid and umbilical cord. **IXIUM TWIN** is a medical device for aid during surgery of the anterior segment of the eye whose rheological characteristics are perfectly suitable for all the various phases in cataract surgery procedures.

IXIUM TWIN :

- maintains the endo-ocular space of the anterior segment of the eye and maintains tissue integrity
- has excellent rheological properties facilitating capsulorrhesis and intraocular lens insertion
- enables excellent visibility of the operative space
- is easy to remove from the anterior chamber
- does not interfere with the cicatrization process
- is non-antigenic and well tolerated by the human eye.

COMPOSITION :

Ingredients	Percentage formula	Function
Sodium hyaluronate (phase 1)	2000 mg	high viscosity
Sodium hyaluronate (phase 2)	1400 mg	low viscosity
Sodium chloride	900 mg	isotonicity
Water for injections	qsp 100 mL	dissolution

INDICATIONS :

IXIUM TWIN is indicated as a surgical aid (medical device) in surgical procedures for cataract, including capsulorrhesis, phacoemulsification, extraction of the crystalline lens and intraocular lens insertion.

IXIUM TWIN maintains the depth of the anterior chamber throughout the duration of surgery enabling enhanced precision of operative procedures with no risk of trauma to the endothelium or other endo-ocular tissues.

SPECIAL PRECAUTIONS FOR USE :

The following precautions for use are recommended during anterior segment surgery :

- check the integrity of the individual sterility protection of the product prior to use ;
- the canula and syringe are single-use disposable for intraocular injection only ;
- the quantity injected into the anterior chamber must be proportional to the volume of the aqueous humor and the anatomical structure to be protected ;
- eliminate all the product by irrigation and/or aspiration at the end of the procedure. A mechanical blockage of drainage at trabecular level may occur, giving rise to a transient increase in post-operative intraocular pressure ;
- the product is to be administered with caution and under special monitoring in patients presenting with pre-existing glaucoma and in the event of glaucoma surgery and/or surgery combined with crystalline lens extraction. In the event of a supra-normal post-operative intraocular pressure, prescribe appropriate treatment ;
- all post-operative inflammatory reactions (iritis, hypopyon, uveitis) and corneal decompensations of the edema type are inherent in surgery of the anterior segment of the eye and no causal relationship with the product has been demonstrated.

INCOMPATIBILITIES :

Do not use quaternary ammonium compounds (benzalkonium chloride) with **IXIUM TWIN**. Sodium hyaluronate precipitates in the presence of quaternary ammonium compounds.

CONTRA-INDICATIONS :

There is no contra-indication to use of **IXIUM TWIN** in compliance with the user package leaflet.

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