
|           |                               |           |                                 |
|-----------|-------------------------------|-----------|---------------------------------|
| <b>EN</b> | <b>INSTRUCTIONAL FIGURES</b>  | <b>FR</b> | <b>CHIFFRES INSTRUCTIONNELS</b> |
| <b>DE</b> | <b>ANLEITUNGSABBILDUNGEN</b>  | <b>ES</b> | <b>CIFRAS INSTRUCCIONALES</b>   |
| <b>PT</b> | <b>FIGURAS INSTRUCCIONAIS</b> | <b>IT</b> | <b>FIGURE DIDATTICHE</b>        |
| <b>RU</b> | <b>ИНСТРУКЦИОННЫЕ РИСУНКИ</b> | <b>EL</b> | <b>ΕΙΚΟΝΕΣ ΔΙΔΑΚΤΙΚΗ</b>        |

|           |                |
|-----------|----------------|
| <b>EN</b> | <b>ENGLISH</b> |
|-----------|----------------|

#### Preloaded Hydrophobic Intraocular Lens

##### DESCRIPTION

Orizon intraocular lens (IOL) is designed for cataract surgery as an anterior surface aspheric optic device to replace the crystalline lens and to reduce the spherical aberration in the human eye. This foldable IOL consists of hydrophobic acrylic manufactured with bonded UV-absorber. It is sterile, one-piece and suitable for posterior chamber implantation. It has a biconvex optic with supporting haptics in *Modified L-Loop* configuration.

##### PHYSICAL PROPERTIES

Material: Hydrophobic Acrylate-Methacrylate Copolymer  
Water content: <1% A-Constant: 118.5 YAG Laser: Compatible  
Optical Design: Aspheric, Biconvex, Square-edge Haptic Angle: 0°  
Sterilization Type: Ethylene Oxide  
Available Diopters: 0 to 10: 1D steps / 10.5 to 30: 0.5 D steps

Orizon IOLs are supplied preloaded within a disposable sterile injector system. The injector system is made of high-quality materials and has been designed exclusively for foldable intraocular lens (IOL) implantation. The injector system allows the ophthalmic surgeon to safely inject Orizon preloaded intraocular lenses into the eye through <2.2 mm incisions.

##### INDICATIONS

Orizon foldable posterior chamber IOLs are indicated for primary implantation for the visual correction of aphakia in adult patients within whom a cataractous lens has been removed. These lenses are intended for placement in the capsular bag.

##### CONTRAINDICATIONS

To date there are no absolute contraindications to intraocular lenses' implantation. Relative contraindications include some forms of chronic acute uveitis and retinal diseases in which the implant may interfere with retinal surgery

##### COMPLICATIONS

Cataract surgery, with or without lens implantation might be associated with: corneal inflammation; hemorrhage; intraocular pressure elevation; post-operative infection; retinal breaks and detachment; cystoid macular edema; corneal edema; posterior capsule opacity.

Complications related with intraocular lens implantation: Capsular rupture; vitreous loss; lens decentration and luxation; wrong calculation of IOL power; damage of IOL during implantation.

As with any surgical procedure, there is risk involved. Potential complications accompanying cataract or implant surgery are not limited to the abovementioned complications.

##### WARNINGS

- Do not use the intraocular lens in case the sterile bag was opened or damaged or in any case of doubt.
- Do not use the intraocular lens after expiration date.
- Do not resterilize this intraocular lens by any method.
- Do not reuse the lens. Single use only. The lens must be used once only for a single patient. Reuse and/or resterilization may compromise device performance and could cause serious harm to patient's health and safety.
- Do not store at temperatures under 5°C and above 35°C. The manufacturer recommends the storage and usage of lenses at room temperature. Do not expose the lens to sunlight.
- Do not soak the lens and injector system in solutions other than balanced saline solutions or equivalents.
- Handle the lens carefully to avoid damage to lens surfaces or haptics.
- Do not attempt to reshape haptics in any way.
- After implantation the lens position must be reversed-S.

##### CALCULATION OF LENS POWER

Preoperative calculation of required lens power for these posterior chamber intraocular lenses should be determined by the surgeon's experience, preference, and intended lens placement. Lens power calculation methods are described in the references found at the end of the text.

##### DIRECTIONS FOR USE

The injector system with preloaded lens is used to fold and implant Orizon IOLs in the capsular bag. **The appropriate surgical technique is the responsibility of the surgeon. He or she must assess the suitability of the given procedure based on his or her education and experience.**

Follow the steps as described below to ensure safe handling.

Prior to implanting, check the lens package for proper lens model, dioptic power, and expiration date.

- Open blister containing sterile injector and the lens and place it in the sterile field. The lens is preloaded on the lens holder section *[see Fig.1]*.
- Hold injector around its midsection without exerting any pressure on the plunger. Rinse the injector nozzle with BSS. Next, fill injector nozzle with viscoelastic material either from above or through the nozzle tip, ensuring correct lubrication of inner nozzle surface. Lubricate both sides of the folding element with viscoelastic *[see Fig.2]*.
- To load the lens from the lens holder in the injector, move the sliders to fold the haptics in the proper position to the end in the direction of the arrow. Check that the haptics have been folded correctly *[see Fig.3]*. *Caution:* Do not continue implantation until haptics are correctly folded onto the lens.
- Gently pull folding element into lens holder until a *"CLICK"* is heard *[see Fig.4]*. This signals that IOL injector system is securely 'loaded' and ready for injection.
- Once loaded, the IOL should be injected immediately, as the viscoelastic may lose lubricity if allowed to stand too long while exposed to air.
- Gently advance the plunger in a slow, even, and controlled manner, while ensuring the IOL haptics are correctly folded onto the IOL optic *[see Fig.5]*. If haptics are caught between plunger and injector nozzle, slacken pressure on plunger until IOL is retracted. *Caution:* Do not continue implantation until haptics are correctly folded onto the IOL.
- To proceed inject IOL by guiding injector tip through the incision *[see Fig.6]*. Stop pressing plunger when IOL exits injector. Do not let the injector soft tip exit the nozzle. Withdraw injector from eye.
- Thoroughly remove viscoelastic material from eye and IOL with standard irrigation and aspiration techniques.
- Disposal: After use, the injector system must be disposed in accordance with your clinical guidance for disposal of single use surgical instruments.

##### PACKAGING

Orizon lenses are packaged in a sterile blister. The packaging sterility is guaranteed until expiry date if blister is not damaged or opened.

##### NOTES:

The lenses are carefully checked and inspected by the manufacturer to assure a high-quality product. Before implantation, study them carefully under surgical operation microscope. If a defect or deformation is noted or suspected, the lens should be returned to the manufacturer. The patient card included in the package is to be completed and given to the patient, together with instructions to keep the card as a permanent record to be shown to any eye care practitioner that the patient may consult in the future. Implantation has to be performed by a trained ophthalmic surgeon to adult patients with aphakia. Orizon posterior chamber lenses are designed for implantation in the capsular bag by injector method through a 2.20 mm incision.

##### SUGGESTED A-CONSTANT

The suggested A-constant indicated on the outer label is presented as a guideline and is a starting point for implant power calculations. It is recommended that you develop your own constant appropriate for you based on clinical experience with the particular lens model, surgical techniques, measuring equipment, and postoperative results. If additional information on lens power calculation is needed, please contact the manufacturer using data given below.

##### DISCLAIMER

The manufacturer will not be liable for any injury suffered to a patient as a result of: any implantation method or technique used by a surgeon to implant the lens; any prescription selection and use of the lens for any individual patient or patient's condition. The manufacturer makes no expressed or implied warranties in connection with the sale of this intraocular lens. The manufacturer is liable only for the design and production of the intraocular lens and not for any incidental, consequential, indirect or exemplary damages of any kind, directly or indirectly arising from the purchase or use of this product.

|           |                |
|-----------|----------------|
| <b>DE</b> | <b>DEUTSCH</b> |
|-----------|----------------|

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