MONOBLUE NafX

CE 0459

Mode d'emploi

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MONOBLUE NafX

Composition - Properties

Purified trypan blue 0.15%, isotonic buffered vehicle
- monopotassic phosphate
- disodic phosphate
- sodium chloride
- mannitol
- deuterium oxide
- water for injection

pH 7.25±0.25, osmolarity 300±50 mOsm/Kg H2O, density adapted.

Description

MONOBLUE NAFX is a purified trypan blue isotonic sterile and apyrogenic solution presented in single dose syringe containing 0.75ml of 0.15% solution. The syringe together with a connector is packaged in a pouch. It is externally sterile. Each box contains 5 syringes **MONOBLUE NAFX** and 5 tuberculin syringes 1 ml. A set of patient labels is provided to allow traceability. **MONOBLUE NAFX** does not contain latex, phthalates or preservatives.

Indications

MONOBLUE NAFX is indicated for use as an aid in ophthalmic surgery by staining the epiretinal membranes and internal limiting membrane during the surgical vitrectomy procedures, facilitating removal of the membranes.

MONOBLUE NAFX is intended to be applied directly on areas where membranes could be present, staining any portion of the membrane which comes into contact with the dye. The dye does not penetrate the membrane.

Contraindications - Precautions

MONOBLUE NAFX is well tolerated by the ocular tissues.

The intraocular tissues integrity should be taken in consideration to control the staining of the target tissue and to avoid non expected staining of other tissues.

Avoid subretinal injection.

MONOBLUE NAFX should be used only for the retinal surgery application.

Do not use in pregnant women or with presumption of pregnancy. Do not use in nursing women or infants. Do not use in patients with known hypersensitivity. Use **MONOBLUE NAFX** with caution when the patient is implanted with a hydrophilic acrylic intraocular lens, the hydrated material having shown potential reversible absorption of the dye.

MONOBLUE NAFX is contraindicated when a non-hydrated (dry state) hydrophilic acrylic intraocular lens is planned to be inserted into the eye, the dye may be absorbed by the lens and stain it.

Thoroughly rinse the excess product with a saline solution immediately after application.

Single use medical device, do not resterilize.

Multiple use of **MONOBLUE NAFX** syringe will not allow guarantying the sterility of the solution any more.

The integrity of the device is guaranteed if the packaging is not damaged. Do not use if the pouch is compromised.

Do not use after expiration date.

Only use cannula between 20G and 25G.

Storage conditions

Store in the original packaging, at ambient temperature, protected from light, in a clean dry place. Do not refrigerate. Do not freeze.

Instructions for use

Prepare the injection syringe: adapt the connector onto the glass syringe. Then connect the tuberculin syringe to the other end of the connector, slowly transfer the total volume of **MONOBLUE NAFX** into the tuberculin syringe by pushing the glass syringe plunger, disconnect the connector and replace by the injection cannula.

Verify the good mobility of the plunger before use, to avoid an unexpected jet which may lead to a subretinal passage of the solution. Protection of the retina may be achieved with the interposition of an instrument between the retina and the cannula output.

Before injection of **MONOBLUE NAFX** a fluid-air exchange is not required, however it may be performed by filling the entire vitreous cavity with air, to prevent aqueous dilution of **MONOBLUE NAFX**. Carefully apply **MONOBLUE NAFX** using a blunt cannula onto the retinal membrane by gently pushing the syringe plunger. Make sure the cannula does not contact the retina to prevent damage. Remove all excess dye from the vitreous cavity before performing air-fluid exchange to prevent the spreading of the dye. **MONOBLUE NAFX** can also be injected in a BSS filled vitreous cavity: 30 to 60 seconds of application under BSS is usually sufficient to achieve appropriate staining.

Waste disposal :

With household waste or recycle	With hospital waste
Paper box, Instructions For Use leaflet, traceability labels and peelable pouches	Glass syringe, connector and tuberculin syringe

Method of sterilization Autoclave

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