Instructions for Use

Capsular Rings

Morcher GmbH

Kapuzinerweg 12 | 70374 Stuttgart | Germany

www.morcher.com

CE 0483

	[en] ENGLISH INSTRUCTIONS FOR USE CAPSULAR RINGS	4
	[bg] БЪЛГАРСКИ ИНСТРУКЦИИ ЗА УПОТРЕБА КАПСУЛАРНИ ПРЪСТЕНИ	17
	[cs] ČESKÝ NÁVOD K POUŽITÍ KAPSULÁRNÍ PRSTENCE.....	31
	[da] DANSK BRUGSANVISNING KAPSELRINGE	44
	[de] DEUTSCH GEBRAUCHSANWEISUNG KAPSELRINGE	57
	[el] ΕΛΛΗΝΙΚΑ ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ ΠΕΡΙΦΑΚΙΚΟΙ ΔΑΚΤΥΛΙΟΙ	71
	[es] ESPAÑOL INSTRUCCIONES DE USO ANILLOS CAPSULARES	85
	[et] EESTI KEEL KASUTUSJUHEND KAPSLIKUJULISED RÕNGAD	99
	[fr] FRANÇAIS NOTICE ANNEAU CAPSULAIRE	112
	[hr] HRVATSKI UPUTE ZA UPORABU KAPSULARNI PRSTENOVİ	126
	[hu] MAGYAR HASZNÁLATI ÚTMUTATÓ KAPSZULÁRIS GYŰRŰK.....	139
	[is] ÍSLENSKA NOTKUNARLEIÐBEININGAR HÝÐISHRINGIR.....	152
	[it] ITALIANO ISTRUZIONI PER L'USO ANELLI CAPSULARI.....	165

	[lv] LATVIEŠU LIETOŠANAS INSTRUKCIJA KAPSULAS GREDZENI	178
	[nl] NEDERLANDS GEBRUIKSAANWIJZING CAPSULAIRE RINGEN	191
	[pl] POLSKI INSTRUKCJA UŻYTKOWANIA PIERŚCIENIE KAPSUŁOWE	205
	[pt] PORTUGUÊS INSTRUÇÕES DE UTILIZAÇÃO ANÉIS CAPSULARES.....	219
	[ro] ROMÂNĂ INSTRUȚIUNI DE UTILIZARE INELE CAPSULARE	233
	[sk] SLOVENSKY NÁVOD NA POUŽITIE KAPSULÁRNE KRÚŽKY	247
	[sl] SLOVENŠČINA NAVODILA ZA UPORABO KAPSULARNI OBROČKI.....	260
	[sv] SVENSKA BRUKSANVISNING KAPSULÄRA RINGAR	273

MODELS

REF

10C, 13, 13A, 13B, 13S, 14, 14A, 14C, 20C, ACTR10, ACTR11, ACTR12, 1G, 1L, 1S, 2C, 2L, 2S, 6D, 6E, 6F, 9L

DISCLAIMER

Please note: Not all models are available in all countries and from all distributors.

PRODUCT DESCRIPTION

Capsular Rings, hereinafter also referred to as implants, are sterile, non-optical medical devices designed to stabilise the capsular bag in patients undergoing cataract extraction. Capsular Rings are made from one piece of clear or blue polymethylmethacrylate (PMMA). They are available with different diameters and in various forms in order to ensure an optimal fixing of the implants in the eye (normal or myopic eye).

Capsular Rings are available in different modifications:

- Capsular Tension Rings (CTR)
- Cionni Rings for scleral fixation
- Cionni Segments for scleral fixation

DEVICE MATERIALS

The clear Capsular Rings consist of > 99 % PMMA and < 0.5 % UV absorber. The blue Capsular Rings consist of > 99 % PMMA; < 0.5 % UV absorber; and < 0.5 % dye.

SPECIFICATIONS

Model (REF)	Remarks	No. of fixation eyelets	Size open	Size compressed	ACCC
10C	CTR	n.a.	12.1 mm	11.0 mm	n.a.
13	CTR	n.a.	12.3 mm	10.0 mm	n.a.
13A	CTR	n.a.	13.0 mm	11.0 mm	n.a.
13B	CTR	n.a.	14.5 mm	12.0 mm	n.a.
13S	CTR	n.a.	10.5 mm	9.0 mm	n.a.
14	CTR	n.a.	12.3 mm	10.0 mm	n.a.
14A	CTR	n.a.	14.5 mm	12.0 mm	n.a.
14C	CTR	n.a.	13.0 mm	11.0 mm	n.a.
20C	CTR	n.a.	12.0 mm	11.0 mm	n.a.
ACTR10	CTR	n.a.	12.3 mm	10.0 mm	n.a.
ACTR11	CTR	n.a.	13.0 mm	11.0 mm	n.a.
ACTR12	CTR	n.a.	14.5 mm	12.0 mm	n.a.
1G	Cionni Ring	1	12.2 mm	10.0 mm	6.0 mm
1L	Cionni Ring	1	13.0 mm	11.0 mm	7.0 mm
1S	Cionni Ring	1	10.7 mm	8.5 mm	6.0 mm
2C	Cionni Ring	1	13.0 mm	11.0 mm	7.0 mm
2L	Cionni Ring	2	13.0 mm	11.0 mm	7.0 mm
2S	Cionni Ring	2	13.0 mm	11.0 mm	5.0 mm

ACCC = anterior continuous curvilinear capsulorhexis; n.a. = not applicable; no. = number

Model (REF)	Remarks	No. of fixation eyelets	Total diameter	ACCC
6D	Cionni Segment	1	9.5 mm	5.0 mm
6E	Cionni Segment	1	10.0 mm	5.0 mm
6F	Cionni Segment	1	7.5 mm	5.0 mm
9L	Cionni Segment	2	11.3 mm	4.8 – 5.7 mm

ACCC = anterior continuous curvilinear capsulorhexis; no. = number

INTENDED PURPOSE

Capsular Rings are intended for the stabilisation of the crystalline lens capsule in the presence of weak or partially absent zonules in patients undergoing cataract extraction with intraocular lens implantation.

INTENDED USER

Capsular Rings are used by either an ophthalmic surgeon or an ophthalmologist during a state-of-the-art operation.

USE ENVIRONMENT

Capsular Rings are used under aseptic conditions in an operating room or a doctor's practice using a surgical microscope.

TARGET GROUP

Patients with fully developed eyes and weak or partially absent zonules undergoing cataract extraction with IOL implantation. Indications associated with weak or partially absent zonules are listed in the section „indications“.

The implants have not been clinically investigated in pregnant or nursing women. Ophthalmologists/ophthalmic surgeons, considering implantation in such patients, should carefully assess the benefit-risk ratio before implantation.

INDICATIONS

Capsular Rings are used only in connection with a cataract operation with IOL and only in connection with IOLs conform to the international standard Ophthalmic implants – Intraocular lenses – (ISO 11979 series).

Capsular Rings help to stabilize the capsular bag in the case of

- high myopia
- zonulolysis
- pseudoexfoliation
- primary zonular weakness
- defective or missing zonules

ADDITIONAL INDICATIONS ONLY FOR CAPSULAR TENSION RINGS

- Marchesani Syndrome

ADDITIONAL INDICATIONS ONLY FOR CIONNI RINGS

- Secondary zonular weakness (e.g. Marfan's Syndrome)

CONTRAINDICATIONS

The usual contraindications for a cataract operation apply here, as well as the following ones:

- patients, whose eyes are not fully developed, i.e. implants are contraindicated in eyes that are still in their developmental stage
- pre-operative complications of the cataract surgery (vitreous humour prolapse, bleeding)
- expected post-operative problems caused by the cataract surgery
- severe microphthalmia or macrophthalmia
- flat front section, e.g. microphthalmia or certain types of chronic angle closure glaucoma; flat anterior chamber
- chronic uveitis
- ICCE (intracapsular cataract extraction) technology
- progressive ophthalmic diseases (proliferative diabetic retinopathy, uncontrolled glaucoma)

Certain operative complications could also lead to a contraindication against the use of the implant, e.g.:

- rupture of the capsule of the lens with or without vitreous prolapsed
- damage to the capsular bag
- prolonged bleeding or other visual impairments
- cortical remnants
- vitreous problems
- posterior capsular opacification (at least in conventional capsular rings without sharp edge)
- IOL decentered
- elevated intraocular pressure (IOP)
- retinal detachment
- deposits on IOL
- iritis
- posterior synechiae

INTENDED CLINICAL BENEFIT

The use of Capsular Rings will lead to the benefit of a stabilised crystalline lens capsule in the presence of weak or partially absent zonules in patients undergoing cataract extraction. This improved anatomical condition of the capsular bag facilitates the technical success of the cataract surgery and contributes to the stabilisation of the IOL and to the improvement of visual acuity.

PERFORMANCE CHARACTERISTICS

Depending on the clinical study, technical success rates of cataract surgery of 92.7 – 100 % and IOL centration rates of 80 – 100 % were reported in the literature. IOL centration was achieved in 94 % and 96.2 % in clinical studies performed by Morcher. For details see section “Summary of safety and clinical performance”.

LIFETIME

The lifetime of the Capsular Rings is 10 years. Without any complications or adverse events, the products may stay in the patient for his lifetime.

COMPLICATIONS OF THE IMPLANTATION (SPECIFIC AND NON-SPECIFIC)

- angle closure
- suture breakage (for sclera-fixated Capsular Rings)
- typical complications of cataract surgery including
 - anterior chamber inflammation (occurring in nearly half of patients)
 - macular oedema
 - corneal oedema
 - corneal dystrophy
 - pupillary block
 - iridocyclitis, vitritis, uveitis
 - precipitates on the intraocular lens

- decentration of the intraocular lens
- (sub-) luxation of the intraocular lens
- infection, endophthalmitis, and panophthalmitis
- hyphema
- elevated IOP
- glaucoma
- vitreous bulge
- vitreous detachment
- retinal detachment
- post-cataract membrane/secondary membrane (after cataract)/posterior capsule opacification (PCO)
- inflammatory pupillary membrane
- iris atrophy
- anterior capsular contraction syndrome
- toxic anterior segment syndrome
- descemet membrane tear and detachment
- anterior synechia
- dysphotopsia

SIDE EFFECTS

The possible side effects of the Capsular Rings are as follows:

- zonular trauma
- impairments caused during cortical cleaning
- possible dislocation of the IOL or the capsular bag when the posterior capsular surface tears

WARNINGS

- Do not use the implant if the sterile packaging is already open or damaged.
- Do not re-sterilize the implant. Re-sterilization can lead to damage to the material and loss of functionality of the implant.
- Do not re-use the implant. It is intended for single use. Re-use can lead to reduced effectiveness or functionality of the implant and/or loss of sight.
- Rinse the implant only with sterile rinsing solutions such as sterile Ringer's solution or sterile BSS® (balanced salt solution).
- Observe the accompanying documents and the warning instructions before using the device.
- Consult the Instructions for Use before using the device.
- Eyes with pseudoexfoliation syndrome and flat anterior chamber run a major risk of zonular defects or other intraoperative complications.
- Do not use a damaged implant. An implant that has been damaged during implantation should not be allowed to remain in the eye.
- If a posterior capsulorhexis (PCCC) is intended by the surgeon, the Capsular Ring must be used after the PCCC. Otherwise there is a risk that the tension of the Capsular Ring may cause the capsulorhexis to lacerate uncontrollably.
- Do not use the implant if the expiration date has passed.
- Do not use the implant if it may have been contaminated with unsterile material (e.g. if it falls on the floor). Dispose of the potentially contaminated implant and use a new product.

PRECAUTIONARY MEASURES

QUALIFICATION AND TRAINING FOR USERS

- The surgeon performing the implantation should carefully read all the information material provided by the manufacturer for the proper use and insertion of the implant before carrying out the implantation.
- The implant may only be used by a doctor experienced in the implantation of intraocular implants, or under his/her supervision. The use of the Capsular Rings requires a detailed understanding of the technical fundamental principles, clinical applications and the risks of the devices.

NOTE WHEN USING LASERS

The laser beam must be focused directly on the target point behind the implant. If a laser beam is focused on the implant, it can result in damage to the implant.

STERILISATION AND PACKAGING

- The implant is sterilised by gamma irradiation.
- The implant is supplied dry and sterile.
- The primary container is sealed, germproof, and packed into a sealed and germproof transparent pouch (double sterile barrier system).
- Sterility is only guaranteed as long as the transparent pouch and the primary container are neither opened nor damaged.
- The expiration date is printed on the folded box and the primary container.

STORAGE CONDITIONS / TRANSPORT CONDITIONS

- Store flat at room temperature + 15 °C up to + 25 °C.
- Transport at - 18 °C up to + 60 °C.
- Protect from direct sunlight and extreme temperatures.
- Protect from moisture/water and store in a dry environment.

PATIENT INFORMATION

- Each patient shall be fully informed about the implant, especially about potential risks (including any undesirable side effects, warnings, precautions, and contraindications) and benefits, before implantation.
- The surgeon performing the implantation shall inform the patient of all side effects and risks, particularly those that have not already been pointed out in the initial consultation.
- The patient shall be instructed to inform the doctor about all side effects, which he/she may observe after the implantation, especially in case of malfunction of the device or changes in its performance.

- The patient shall be informed that there are no specific environmental conditions known that can affect the Capsular Rings. Medical procedures or examinations like laser treatment, magnetic resonance imaging (MRI) or X ray scan do not negatively impact the performance of the Capsular Rings and do not pose any risk to the Capsular Rings.
- The patient shall be informed that there are no known circumstances, in which the Capsular Rings interfere with certain diagnostic tests, evaluations, or therapeutic treatments or other procedures, or interfere with the equipment used in these procedures.
- Written information about the devices intended for the patient is available at <https://patient-information.morcher.com>

IMPLANT CARD

The relevant details and one label shall be attached onto the enclosed implant card. The date of the implantation shall be written on this card and the eye that was operated on shall be crossed off. Instructions for filling out the implant card are available at <https://ifu.morcher.com>.

The completed card shall be given to the patient for safekeeping to ensure traceability of the implant. It should be presented to every eye specialist whom the patient consults.

NOTES ON HANDLING / IMPLANTATION

1. Check all labels with regard to type, specific details and expiration date before opening the transparent pouch.
2. Only pick up the implant with smooth-tipped forceps and handle it carefully to avoid damage.
3. Under aseptic conditions, remove the implant as follows:
 - open the transparent pouch at the points indicated and remove the primary container
 - carefully pull off the seal from the primary container
 - carefully remove the implant; do not pull on the implant
4. Check the implant before use and only use an undamaged implant.
5. It is recommended to rinse the implant with sterile Ringer's solution or sterile BSS® prior to the implantation.
6. Implantation:
 - Before implantation, fill the capsular bag with sufficient viscoelastic to prevent the anterior and posterior capsules from sticking to one another and thus facilitate insertion of the ring into the capsular bag.
 - When using forceps to implant the Capsular Ring, dial the ring through the incision and the capsulorhexis into the capsular bag with a holding forceps.

When using an injector to implant the Capsular Ring, we recommend the Geuder® injector G 32960. Follow the injector manufacturer's instructions for use to avoid problems during implantation and damage to the implant. Capsular Tension Rings should be pulled into the injector completely; the fixation loops of Cionni Rings must remain outside the injector's cannula to avoid damage to the implant.

- We recommend inserting the ring with an acute angle in a tangential direction to avoid radial pressure from the leading eyelet on the equator. The Capsular Rings are designed in such a manner that they follow the natural curve of the capsular bag.
- If a posterior capsulorhexis (PCCC) is intended by the surgeon, the Capsular Ring must be used after the PCCC. Otherwise there is a risk that the tension of the Capsular Ring may cause the capsulorhexis to lacerate uncontrollably.

SELECTION OF A PROPER CAPSULAR RING

CTR are usually used in cases of mild zonular instability (i.e., zonular dehiscence < 4 clock hours); in cases of progressive or advanced zonulopathy Capsular Rings for scleral fixation should be used.

Larger Capsular Rings may be used in larger capsular bags; axial length and corneal diameter correlate with the size of the bag and may be used for Capsular Ring size selection. However the use of a larger ring may be chosen since overlap of the end terminals provides a complete circumferential support, although it may be more challenging to insert.

TIMING OF CAPSULAR RING INSERTION

A Capsular Ring may be inserted at any stage of the cataract procedure. We recommend to insert Cionni Rings after nuclear and cortical removal.

INSTRUCTION FOR IMMEDIATE IMPLANTATION OF AN INTRAOCULAR LENS

The cradle of the implanted open loop lens or the axis of the IOL must be set at a 90° angle to the axis of the capsular ring. It is an advantage to have the cradle of the lens positioned posterior to the ring in the fornix of the capsular bag. This can be achieved by applying pressure on the lens in the region of the haptic or on the cradle itself.

ADDITIONAL INSTRUCTIONS FOR CIONNI RING

The size of the anterior rhexis depends on the implant used.

When inserting a Cionni Ring, ensure that the leading fixation loop (the curved end of the ring) with the eyelet is located outside of the capsular bag.

Position the fixation eyelets of the Cionni Ring at the position of the zonular defect.

Once the Cionni Ring is in the correct position, it should be sutured to the sclera using the fixation eyelet. A double-armed 9-0 prolene suture should be used for fixation of the implant.

REMOVAL OF THE RING

The Capsular Ring may stay in the eye for a lifetime. However, it is sometimes difficult to implant the ring satisfactorily into the bag on the first try and removal may be necessary. If it becomes necessary to remove the ring from the capsule at any time, use of a hook to lift one end of the ring out of the capsulorhexis is recommended. Then radial the ring out of the bag. In some situations it may be desirable to cut the ring into two pieces with vannas scissors and to dial each piece separately out of the bag. Another option is the use of a CTR injector to withdraw the CTR in one piece through the initial phaco incision.

Note: When removing the ring the sutures must be cut and removed as well.

USE WITH OTHER DEVICES AND / OR GENERAL PURPOSE EQUIPMENT

The Capsular Rings are to be used in connection with IOLs conform to the international standard Ophthalmic implants – Intraocular lenses – (ISO 11979 series).

There are no restrictions known regarding the combination of the Capsular Rings with ophthalmic viscosurgical devices (OVD, viscoelastics), intraocular physiological rinsing solutions, or standard surgical equipment/instruments. When using an injector to implant the Capsular Ring, we recommend the injector G - 32960 made by Geuder® AG, Heidelberg, Germany.

For scleral fixation of Cionni Rings and Ring segments double-armed 9-0 prolene sutures should be used as these have proven long-term stable.

DISPOSAL

- The packaging and its contents shall be disposed of in accordance with the current local regulations.
- There are no special disposal regulations applicable for devices, which are uncontaminated.
- Contaminated devices/packaging shall be handled and disposed of in accordance with the provisions of the respective clinic.

EXCLUSION OF LIABILITY

Morcher GmbH and its distributors cannot be held liable for any procedures, methods, or surgical techniques used by the surgeon, nor for the selection of the implant for any eye. This is solely the responsibility of the surgeon.

COMPLAINTS

Complaints should be addressed to

Morcher GmbH

Quality Assurance Department
Kapuzinerweg 12
70374 Stuttgart
Germany

Phone + 49 (0) 711 / 95 320 - 0
Fax + 49 (0) 711 / 95 320 - 80
E-Mail complaints@morcher.com

or to your distributor.

REPORTING ON ADVERSE REACTIONS

Any serious incident that occurred in relation to the device should be reported to

Morcher GmbH

Quality Assurance Department
Kapuzinerweg 12
70374 Stuttgart
Germany

Phone + 49 (0) 711 / 95 320 - 0
Fax + 49 (0) 711 / 95 320 - 80
E-Mail complaints@morcher.com

and to the competent authority of the member state in which the user and/or patient is established.
Please also inform your distributor if necessary.

RETURN OF A DAMAGED / NON-STERILE DEVICE

If an implant is returned, the serial number, batch code and reason for its return must be provided. The packaging, containing the returned item, must be in accordance to the condition of the implant and, where necessary, must warn of any contamination.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The summary of safety and clinical performance (SSCP) referred to in Article 32 of the Medical Device Regulation (EU) 2017/745 is linked to the Basic UDI-DI number 4048509CR01VM and is available at <https://ifu.morcher.com> or in the European database on medical devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed> (as soon as this EUDAMED functionality is available).

ABBREVIATIONS AND PACKAGING SYMBOLS

ABBREVIATION / SYMBOL	DESCRIPTION
IOL	Intraocular Lens
CTR	Capsular Tension Ring
PC	Posterior Chamber
AC	Anterior Chamber
ACCC	Anterior Continuous Curvilinear Capsulorhexis
PCCC	Posterior Continuous Curvilinear Capsulorhexis
ICCE	Intracapsular Cataract Extraction
ECCE	Extracapsular Cataract Extraction
PCO	Posterior Chamber Opacification
D	Diopter in aqueous humour
ØB	Body diameter
ØT	Total diameter
ØO	Optic diameter
ØI	Inner diameter
R	Radius
Th	Haptic side view
A-CON	A-constant
AC DEPTH/ACD	Anterior Chamber Depth
PMMA	Polymethylmethacrylate
PEMA	Polyethylmethacrylate
UV	Ultraviolet
HEV	High Energy Visible light
C-F-M	Compression forged method
CQ	Clinical Quality
BSS®	Balanced Salt Solution
SN	Serial Number

ABBREVIATION / SYMBOL	DESCRIPTION
UDI	Unique Device Identification
GTIN	Global Trade Item Number
MD	Medical Device
REF	Catalogue number
LOT	Batch code
	Use by: Year - Month
EXP	Use by: Year - Month
	Manufacturer
DE	Country and date of manufacture
	Importer
	Distributor
CE XXXX	CE Mark with 4-digit number of the notified body
STERILE R	Sterilised using irradiations
STERILE	Sterilised using steam
	Double sterile barrier system
	Single sterile barrier system
kGy	Kilogray
	Do not reuse
	Do not resterilize
	Consult instructions for use
 ifu.morcher.com	Consult electronic instructions for use
	Caution, consult accompanying documents
	Keep dry
	Keep away from sunlight
	Temperature limit (The applicable lower and upper temperature limits for storage are indicated on the device packaging. The applicable lower and upper temperature limits for transport are indicated on the outer packaging used for transportation.)

ABBREVIATION / SYMBOL	DESCRIPTION
	Do not use if packaging is damaged
	This side up
	Recyclable packaging
	Date of implantation
	Healthcare institution
	Patient information website
	Patient identification

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